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

OPERATOR'S GUIDE

NORLAND

NORLAND XR-series X-RAY BONE DENSITOMETERS
Version 2.2 Software Upgrade
January 15, 1991

Please place these pages at the front of your existing XR-series Operator's Guide. They include instructions for the use of features new to Version 2.2 software. Note that the last page of this update package replaces an existing page of your Operator's Guide.

The following describes the new and enhanced Version 2.2 XR-series X-Ray Bone Densitometer Software features.

Image Menu

A new Image menu has been added which permits all on-screen images to be smoothed via interpolation or displayed in their raw format. Smoothing is faster on High-Resolution or Small Images (but provides little image improvement); slower on Low-Resolution or Large Images (but provides significant image improvement).

To activate image smoothing, select **Setup** and click on the **Image** command. The choice of smoothed or raw format image is presented. Make your selection and click on the **[Accept Changes]** button.

Improved <Enter> function in dialog entry

<Enter> now acts like <Tab> during patient or instrument data entry operations (anywhere there is at least one user-enterable text field). Prior software caused the dark-bordered response button to be triggered when <Enter> was pressed, usually an acceptance of entered values and a return to the previous screen. <Enter> now acts like the <Tab> key, moving the data entry cursor to the next available data entry field.

Pressing <Shift> and <Enter> simultaneously acts like clicking on the dark-bordered button.

Faster New Patient and Patient Search routines

New patient addition and patient record locating routines are now much faster. For example, on XR-series host computers having more than 1000 patients, a 400% performance improvement can be seen.

Soon after installation of 2.2 software, the patient records should be sorted by selecting **Disk Utilities** and clicking on the **Optimize Patients File** command. This operation can take 30 or more minutes to organize the records of a very large patient file or only a few minutes for smaller files. After being sorted, future patient access and new patient addition operations will be substantially quicker. This operation can sort existing system, volume, or copy disk Patients Files.

Norland recommends that the hard disk patient records be sorted, as described above, once each month or so, to preserve the enhanced operations.

Improved Calibration Procedures

Calibration procedures have been improved. Should in-process errors be discovered, messages appear detailing errors detected. The system will automatically retry any failed sample. The instrument procedure continues on in event of any error scanning the calibration standard (however files from failed calibration are saved under alternate names, for diagnostic purposes.) And Calibration Results screens now show the x-ray source type and serial number.

Print-Screen

A Print-Screen feature has been added to existing print features and operates on almost every screen in XR-series software. It provides a hard copy of the screen contents on the standard XR-series parallel printer.

To use, simply press the **<Print Screen>** key (upper right on keyboard) for an immediate 'screen dump' to the printer - no queue, one copy - it's slow, but it does replicate the screen exactly (in the displayed color or grayscale mode).

An information box appears, asking if the image is to be printed as is or rotated 90° (to fill the page). Click on the desired choice and the printout is initiated. To preserve instrument operations, do not attempt to interrupt this printout.

Service Setup Menu

A Service Setup menu has been added which brings all source setup and scanner information onto one screen. It is here where your installer would enter values. It is not likely that you would use this screen unless your instrument receives certain hardware changes.

Improved Disk Utility Messages

Minor wording changes were made to Disk Utility screens, improving message clarity and operator instruction. To review some of these changes, select **Disk Utilities** and click on any of the listed commands.

Non-U.S. DOS Protection

The XR-series Software Installation procedure now operates properly when a non-U.S. version of the Disk Operating System is installed. Prior versions of the installation procedure would 'hang' because of errors generated by the non-U.S. DOS program.

Other Improvements

Additions have been made to XR-series software to accommodate high-speed scanner procedures in calibration and scanning operations and to accommodate Optical Disk Drive technology.

And, Version 2.2 XR-series software adds more than 25 code changes to enhance operations which aren't visible to the operator and do not have any effect on calculations or measurement results.

Attached is a new page (pages 117-118) for your basic Operator's Guide. After reading the new CV Specification material, please replace your present page with this new page. Thereafter, whenever discussing or using the CV Specifications, keep this information in mind.

The descriptions and narratives found in this Operator's Guide are intended to provide a foundation for use of the Norland XR-Series X-Ray Bone Densitometers using Version 2.2 or later software. The actual screen images, printout images, and depictions of these images seen this guide may not be identical to the installed version of software and are not intended for diagnostic purposes. All specifications are subject to change without notice. All Rights Reserved.

Document Release Approval			
Engineering:	Service:	Sales/Marketing:	QA/RA:
<i>ED. 1/22/91</i>	<i>KE 1/22/91</i>	<i>QuRS 1/22/91</i>	<i>Jef</i> <i>22 Jan. 91</i>

The following table lists the page changes and revision history of this Base XR-Series X-Ray Bone Densitometer Operator's Guide, NORLAND Part Number 388M402:

<u>Rev.</u>	<u>Date</u>	<u>ECN No.</u>	<u>Pages Affected:</u>
A	7-12-90	6320	Software Version 2.1 Production Release, pages i-vi, 1-134.
B	1-16-91	6475	Software Version 2.2 Production Release, this page and pages 117-118.

XR-26 and XR-26 Mark II X-Ray Bone Densitometer

DETAILED SPECIFICATIONS

SYSTEM

Description	Dual Energy X-Ray Absorptiometer
Installation and Training	Complete on-site installation and comprehensive training included in purchase price.
Warranty and Service	One Year parts and labor from date of installation. Extended warranty contracts available.
Operating Temp.	15-32°C (60-90°F).
Humidity	Up to 80% non-condensing.
Calibration	Automatic with supplied Calibration Standard.

A-P SPINE MEASUREMENT

Scan Site	Lumbar Spine (L2-L3-L4) - Anterior-Posterior View.
In-Vivo Precision ¹	BMD: 1.0% Total L2-L4; 2.0% Individual Vertebra. BMC: 1.5% Total L2-L4; 2.0% Individual Vertebra. Area: 1.0% Total L2-L4; 1.0% Individual Vertebra.
Accuracy ²	BMD, BMC, Area: Within 1.0% of industry standard.
Entry Dose	4.4 millirems maximum, less for body thickness less than 18cm.
Spatial Resolution	1.5mm x 1.5mm: Line spacing x point resolution (pixel size).
Scan Speed	60mm/sec.
Scan Time	5 minutes typical.

HIP MEASUREMENT

Scan Site	Femoral neck, Ward's Triangle, and greater trochanter.
In-Vivo Precision ³	BMD: 2.5% neck; 2.0% trochanter; 4.5% Ward's Triangle. BMC: 2.0% neck; 4.5% trochanter; 4.5% Ward's Triangle. Area: 1.5% neck; 4.5% trochanter.
Accuracy ⁴	Within 1.0% of industry standard.
Entry Dose	8.8 millirems maximum, less for body thickness less than 22cm.
Spatial Resolution	1.0mm x 1.0mm: Line spacing x point resolution (pixel size).
Scan Speed	45mm/sec.
Scan Time	6 minutes typical.

SCANNER

Active Scan Area	198cm x 65cm (78" x 25.5").
Drive System	Stepper motors, timing belt drives, scanning in any direction.
X-Ray Source	Stationary anode X-Ray tube; 0.9mm focal spot size; 100 kV constant potential, 1mA constant anode current; Samarium filter (K-edge = 46.8 keV); Minimum filtration is 3mm Aluminum equivalent.
X-ray Detector	Two NaI scintillation detectors in pulse counting mode.
Positioning Aid	Less than 1mw Class II red helium-neon laser.
Electrical	100/120/220/240 VAC \pm 10%, 50/60 Hz, 600 VA maximum.
Physical Spec.	Table: 239cm (94")L x 84cm (33")W x 66cm (26")H. System: 239cm (94")L x 122cm (48")W x 131cm (51.5")H. Tabletop-to-arm distance: 38cm (15"). Weight: 368kg (810 lbs).

STANDARD CONTROL/ANALYSIS WORKSTATION

Computer	IBM PS/2 Model 55SX with minimum 1 MByte RAM memory, minimum 30M Byte hard disk, 1.44 MByte floppy disk, math co-processor, DOS 3.3 or later, Mouse pointing device.
Display	Zenith ZCM 1490 Flat Screen Color Display or equivalent with VGA high resolution graphics.
Printer	Fujitsu DX2300 or equivalent color printer provides Results Reports in color and/or grayscale.
Software	Version 2.1 or later, resident on hard disk; updates supplied on floppy diskettes. With: Automatic bone edge determination; automatic and manual Region-of-Interest identification; system self-diagnostics; patient bone density Trending.
Physical Specs.	Workstation Size: 121cm (48") L x 76cm (30") W x 71 cm (28")H. Weight: 68 kg (150 lbs).

NOTES


All specifications subject to change without notice

- 1 Based on 4 subjects measured 4.5 times (avg.) over 24 days per procedures outlined in the Operator's Guide.
- 2 Based on scans of an industry-standard Anthropomorphic Spine Phantom.
- 3 Based on 9 subjects measured 4.2 times (avg.) over 24 days per procedures outlined in the Operator's Guide.
- 4 Based on the same analysis algorithms and calibration used for spine measurement. Therefore, traceable to the industry-standard Anthropomorphic Spine Phantom.


900630

RECOMMENDED XR-26 SYSTEM LAYOUT

a	Workstation Height	71cm (28")
b	Scanning Platform Height	66cm (26")
c	Maximum Patient Clearance	41cm (16.3")
d	Overall Scanner Height	131cm (51.5")
e	Minimum Wall Clearance	5cm (2")
f	Workstation Width	76cm (30")
g	Minimum Door Width	89cm (35")
h	Recommended Room Width	366cm (144")
i	Workstation Length	122cm (48")
j	Overall Scanner Width	122cm (48")
k	Minimum Wall Clearance	2cm (1")
l	Overall Scanner Length	239cm (94")
m	Recommended Room Length	366cm (144")

Recommended Working Area - Front View

Recommended Working Area - Top View

1. INTRODUCTION TO THE SYSTEM

The Norland XR-26 Series X-Ray Bone Densitometers use the technique known as Dual Energy X-ray Absorptiometry (DEXA) to estimate skeletal bone mineral *in vivo*. These instruments may be used to make measurements of the lumbar spine, proximal femur (hip), and portions or all of the entire body. They may also be used in research protocols to quantify small subject bone mineral content. The XR-26 Bone Densitometers cannot detect existence or absence of bone disease or bone disorders. They may only aid in estimating changes in bone mineral density (BMD). The clinician alone is responsible for performing diagnoses.

Performance of this device is subject to the Detailed Specifications provided at the beginning of this Operator's Guide. Proper operation of the XR-26 and interpretation of its measurements requires a thorough understanding of the factors affecting precision and accuracy. Such information is provided in Section 5.6, "Precision and Accuracy."

The XR-26 Bone Densitometer consists of two major components: a Scanner Unit capable of positioning the x-ray beam over any part of the patient's body; and a Control/Analysis Computer System, which allows the operator to control the scanning process, analyze and store the resulting data, and print reports of the results.

The Scanner is a specially-constructed table which accommodates the patient in a supine position. The easily removable patient pad is cushioned and upholstered with a heavy-duty fabric. Fixtures are provided for special positioning of the patient during spine, hip, and other procedures.

The Scanner Arm incorporates a touch-sensitive Scanner Control Panel. The operator uses the keypads on this panel to position the scanner arm directly over certain anatomical landmarks to define the desired

scan area. A low-power laser is located in the scanner arm. It indicates exactly the position of the photon beam by projecting a spot of red laser light directly on the patient surface.

The radiation source is an x-ray tube with heavy K-edge filtering. The x-ray tube is mounted within a lead-shielded chamber which has an electrically-operated shutter across the beam exit path. The x-ray source is specially designed for densitometry and has a low operating temperature and long life.

The Model XR-26 incorporates a member of the IBM Personal Computer product line¹. This Operator's Guide describes equipment operation using an IBM-PS/2 Model 55SX computer. Although another member of the IBM-PS/2 family may be supplied with the system, all instrument operational details remain very similar to the commands described herein. A standard IBM 'mouse' pointing device provides one-handed program operation. A color CRT display with full graphics capability provides detailed and informative images. A printer, also with graphics capability, provides hard copy printouts of bone measurement results.

Because of the special needs of the XR-26 programs, the computer may not be able to be used to run other software when it is not being used as a bone densitometer. Please check with your local Norland representative before altering the CONFIG.SYS or AUTOEXEC.BAT programs.

NOTE: Throughout this guide, the term "XR-26" refers to the Norland Model XR-26 or Model XR-26 Mark II X-Ray Bone Densitometer. And, although a female patient is characterized herein when describing instrument operations, the procedures remain the same for a male patient.

¹ IBM, PS/2, PC-DOS are trademarks or registered trademarks of International Business Machines, Incorporated.

1.1 CONTRAINDICATIONS

As with any procedure involving radiation, care must be taken to avoid exposing a developing fetus. Since bone density measurements are not critical to short term patient care, these XR-26 procedures are contraindicated for pregnant subjects.

Other contraindications are based on the probability that measurement results could be compromised, such as:

- Subjects who have recently (within seven half-lives) undergone radionuclide uptake procedures. Residual emission may be misinterpreted by the XR-26 detectors as energy generated by the XR-26 X-ray source.
- Subjects who have recently (within one week) ingested or been injected with a radiopaque substance, such as barium, which could result in increased soft tissue photon absorption. Even calcium supplements may be detected in a spine scan.
- Subjects known to have prosthetic devices, implants, or other sub-dermal metallic objects within the scan area, e.g., hip prostheses, which could result in extreme density measurements.

Screening of patients for any of these conditions should become a normal part of pre-scan procedures.

1.2 RADIATION SAFETY

The x-ray source used in the XR-26 instrument is completely self contained and requires no operator maintenance.

CAUTION: Service of the x-ray source must be performed only by authorized Norland Service Representatives.

1.2.1 SAFETY FEATURES

The XR-26 is designed to be safe for both patient and operator:

- The exit beam is highly filtered and highly collimated so as to minimize patient dose.
- The x-ray beam is present only during the calibration process and during scans. Whenever the system is not actually gathering useful data, the computer closes the x-ray beam electrical shutter. Source current and voltage are turned off upon initiation of analysis or when a scan is cancelled.
- Pressing the readily-accessible **EMERGENCY OFF** key results in immediate shutdown of x-ray power, laser power, and scanning operations.

1.2.2 RADIATION DOSE INFORMATION

The Radiation Source

The radiation produced by this instrument, although highly localized, is quite intense in the primary beam path, and should be regarded with due respect.

The source of radiation is an x-ray tube operated at 100 kilovolts, constant potential and 1.0 milliamps current. The tube is located in a shielded enclosure which reduces external radiation to a safe level everywhere except at the beam exit port. The beam port is fitted with an electrically operated shutter which cuts off the beam at all times except when a scan or a calibration is in progress.

Patient Dose

The radiation dose to the patient is dependent on the type of scan procedure and the body thickness of the patient. This latter effect is because the XR-26 uses variable filtration to assure a transmitted count rate for each patient which is consistent with the precision specification. However, the operator does not select the filters.

In actual operation, a test is performed on each patient just before the scan begins. In effect, the test measures the thickness of the patient, as follows: the XR-26 inserts maximum filtration (attenuation) and measures the count rate; it selects a filter appropriate for patient thickness which will provide the proper count rate; that filter is used for the subsequent scan procedure.

Note that even if this variable filtration is at its minimum, the fixed filtration remains in place, so that the aluminum equivalent specification for the instrument remains valid (see Detailed Specifications).

Table I lists typical entrance skin dosage for each scan type based on the listed body thickness. A reference to the scanning parameters used by the XR-26 for each series of measurements is given. Measurements were made with a Victoreen Model 450P Ion Chamber.

Procedure - (Resolution @ Speed)	Body Thickness (cm)	Entrance Dose (mr)
<u>A-P Spine</u> - 1.5mm x 1.5mm @ 60mm/sec	0 to 9	less than 0.15
	9 to 12	0.3
	12 to 15	0.5
	15 to 18	1.1
	18 and over	2.2
<u>Hip Scout</u> - 1.0mm x 2.0mm @ 80mm/sec	0 to 10	less than 0.04
	10 to 13	0.06
	13 to 16	0.1
	16 to 19	0.2
	19 to 22	0.4
	22 and over	0.8
<u>Hip Measurement</u> - 1.0mm x 1.0mm @ 45mm/sec	0 to 10	less than 0.2
	10 to 13	0.3
	13 to 16	0.5
	16 to 19	0.9
	19 to 22	2.1
	22 and over	4.4

Table I - Entrance Skin Dosage

Operator Dose

The dose to the operator is negligible. Leakage radiation from the x-ray source at the surface of the scanner when in standby mode does not exceed background radiation.

During a scan procedure, the leakage radiation is somewhat higher, but still less than 1 mr/hr. The total radiation level next to the scanner is higher due to scatter from the beam passing through the patient. The level at the edge of the scanner tabletop may be as much as 2.5 mr/hr, the exact amount depending on the thickness of the patient.

At this emission level, an operator sitting adjacent to the scanner will receive a dose of 0.2 mr during a 5 minute scan procedure. By moving a distance of 1 meter from the scanner, this dose could be cut by a factor of at least 4. At that distance, the operator could perform 1000 scans and receive a total dose of 50 mr, which is less than the average annual background dose of 80 mr/yr.

1.3 THE LASER POSITIONING SYSTEM

A low-power (less than 1 milliwatt) Class II helium-neon laser is incorporated into the instrument. The laser is turned ON by pressing the LASER ON/OFF key on the Scanner Arm Control Panel. Each press of this key will turn the laser on or off. When the laser beam is activated, the LASER ON indicator is illuminated. The operator is directed to use the laser at appropriate times during scan procedures.

The laser beam is projected along the same path that the x-ray photon beam follows. Where the patient intersects the laser beam, a small, red spot of light appears. By using the spot of laser light as a guide, the scanner arm may be easily positioned over a selected anatomical location on the patient. That location may then be acknowledged by the control/analysis unit for definition of a soft tissue point or an extent of the scan area.

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength.
DO NOT STARE INTO BEAM.

The scanner arm incorporates a manual laser shutter. When closed, the shutter completely blocks the exit of the laser beam onto the scanner table. Certain XR-26 instruments use a Laser Manual Shutoff paddle, which is elevated to horizontal position to open the laser shutter, and depressed to its down position to close the shutter. The Laser Manual Shutoff is in the laser beam path, not in the x-ray beam path. Therefore, it makes no difference in the measurement values whether the Laser Manual Shutoff is open or closed.

Other XR-26 instruments use a rotating laser shutter. To open the laser shutter, turn the Laser Manual Shutoff knob fully clockwise. To close the shutter, turn the Laser Manual Shutoff knob fully counterclockwise. This shutter must always be fully open when performing measurements.

1.4 OPERATOR INTERFACE

Version 2.0 and later XR-26 software is devised within the pictorial, visually-oriented structure of the GEM®/3² operating environment. This 'mouse'-based³ application is designed to take full advantage of the excellent graphics and information-handling capabilities of GEM®/3. This operating environment is installed in your XR-26 computer and configured for its hardware.

Because of the special needs of XR-26, it may not be possible to run other DOS applications on the system when it's not being used for densitometry. If you desire to use the host computer for other programs, please discuss the matter with your local Norland Representative.

Object-Oriented Paradigm

The central philosophy used to guide the design of the Version 2.0 Operator Interface is that of a patient's visit to a doctor's office. In essence, when a patient first arrives at the office, she announces herself to the receptionist. The receptionist confirms the appointment and if she is a new patient, a new patient file is created for her. Otherwise, the patient's current files are obtained. The patient then moves to the examination room. The ordinary routine of the consultation and diagnostic process takes place, and ends with recommendations and the patient's release.

In a similar fashion, the first task for the XR-26 operator is to acquire the patient's records by selecting the patient from the patient list. If she

² GEM/3 and GEM/3 Desktop trademarks of Digital Research Incorporated.

³ Software application designed to work with a mouse.

is a new patient, a record and patient personal history file is created. If she is a repeat patient, her records are updated. All subsequent patient and scan procedures are performed on the selected patient until another patient is selected. The scan procedure to be performed is requested by selecting the type of scan. Then the scan routine is followed, as prompted by display messages. When the procedure is completed, the scan records are saved. The doctor then makes an analysis of the scan and consults with the patient (or referring physician). When the patient is released, the analysis records are saved and the process is complete. The next patient is selected and attended to in the same fashion.

Although this example office visit procedure may not be used worldwide, it does provide a framework for a basic understanding of the sequence of XR-26 software operations.

1.4.1 MOUSE TECHNIQUES

All program actions are initiated by moving an arrow-shaped pointer (seen at left) around on the screen. The pointer is positioned by moving the mouse on any clean flat surface.

The program pointer, as depicted above, moves in coordination with mouse motion and permits selection of menu commands. To best use the mouse, position the device with the cord pointing away from you.

Figure 1 indicates how to hold and move the mouse. The cord should move freely. Although right-handed operation is shown, the mouse may also be used with the left hand.

The mouse doesn't need a special pad or grid - just a few inches of clear desk space. If you run out of mouse movement space, simply pick up the mouse and set it down in the center of your cleared area. The pointer will not move when the mouse is not in contact with the surface.


Figure 1 - Holding and moving your mouse

It is best to hold the mouse in a relaxed manner, with your thumb and little finger at either side of the mouse. Be relaxed and use a light touch. Rest your wrist naturally on the desk surface. Let your hand pivot from the wrist to move the mouse. The faster the mouse is moved, the faster the pointer moves.

XR-26 mouse operation is automatic: i.e., if the mouse is connected to the computer, it will become the principal command input device. If a mouse is temporarily unavailable, XR-26 provides keyboard equivalents to mouse techniques, as described in Section 1.4.3 below.

The buttons on the mouse are used to tell the XR-26 program that an operation or procedure is desired. If your mouse has more than one button (most do), use only the button on the left. XR-26 does not respond to input from the other buttons.

Table II lists all of the mouse techniques used in XR-26. Examples that use each technique follow.

<u>TECHNIQUE</u>	<u>MOUSE ACTION</u>
<i>Click</i>	Press the button once.
<i>Shift-click</i>	Hold down the Shift key (on the keyboard) as you press the button.
<i>Double-click</i>	Press the button twice rapidly.
<i>Drag</i>	<ol style="list-style-type: none"> 1. Press and hold down the button. 2. Move the mouse. 3. Release the button.

Table II - Mouse Techniques

Several words and phrases are used consistently throughout this Operator's Guide. They have special meaning within XR-26 software and are used only in reference to a command input. These words and phrases are:

- *Click* (Clicking)
- *Drag* (Dragging)
- *Shift-click*
- Select
- *Double-click*
- Continue

Click means press and immediately release the mouse button. You click on an object to use it. The XR-26 program asks that you first select an object, such as Select a Patient, and then do something to that object, either by direct manipulation with the mouse or by choosing a command from a menu.


To de-select a command or remove a drop-down menu from the screen, move the pointer to an open area of the display (not on top of any text or images) and click again. This operation removes the highlight from the menu item. Another click in the open area causes the menu itself to disappear.

Shift-click means hold down the keyboard's Shift key while you click on the desired object. Shift-click is used in XR-26 to select all displayed Region of Interest cursors in order to move them as a group.


Double-click means to rapidly press and release the mouse button twice. This operation takes you to the next logical place in the program without having to select and click on Main Menu bar titles.

For example, double-click on a patient name and the list of scans taken for that patient appears. Double-click on a patient scan, and an analysis results display appears.

Drag means to press and hold down the mouse button while you move the mouse. To move an edge of a cursor, position the pointer within the 'control point' located on the cursor edge. Press the mouse button and hold it. Notice that a 'finger' symbol (similar to that shown at the right) replaces the 'arrow' pointer. Move the mouse and observe that the edge you are dragging is moved as you move the mouse.

To move all region cursors as a group without changing their relative positions or sizes, select the cursor group using the Shift-Click technique. Position the pointer at a convenient location near the Region(s) of Interest, then press and hold the mouse button. Notice that a 'fist' symbol replaces the 'arrow' pointer. Move the mouse and observe that the entire cursor group is dragged along with the fist as you move the mouse. Drag the cursors to their new position and release the mouse button.

Select, in the context of XR-26 operation, means to position the pointer on top of the desired Main Menu bar title. It also means to move the highlight up or down within a drop-down command list.

When selecting a Patient or Scan record, you must first move the highlight to the desired name or record, then click.

Continue, in the context of XR-26 operation, means to move on to the next logical operation or screen. Much of this action is performed

by the program. However, certain operations require extensive and variable operator interaction. Therefore, certain screens contain a button with a **Continue** legend. Click on this button to move on to the next step of the program.

Command Appearance

To progress through each procedure described in this Guide, the operator must interact with messages seen on the display. The following type face conventions are consistently used throughout this Guide.

- When a command is shown in a *shadowed* type face, such as **File** or **Regions**, it signifies a Main Menu bar title or a command line listed in a drop-down Main Menu box.
- Similarly, commands shown in an *outlined* type face, such as **Accept Changes** or **Continue**, signify a button legend or other command which will not appear as a Main Menu bar title or drop-down menu command line.

Command Symbol

The *Close Box* symbol (shown at right) appears at the top left of every screen. Clicking on this symbol causes the program to move automatically to the safe completion of an operation. For example, click on the **Close Box** button when done with defining a region, and the next screen to appear is whichever Results screen was last used.

1.4.2 DROP-DOWN MENUS

XR-26 Version 2.0 software features the use of Drop-down Menus. When you place the pointer on top of one of the Main Menu bar titles (*File, Scan, Analysis*, etc.), a sub-menu 'drops' down. Each line in the drop-down menu is a command.

Some commands are 'dimmed' or 'de-emphasized'. You cannot choose these dimmed commands because: they have no meaning in the context of what operation is in progress at the time; it is already the command in effect; the option is not installed.

The command "Show Baseline" (shown at right) is not available, since the display is already presenting the baseline information. All you can do with baseline now is "Hide" it. You can also see that "16 Level Gray" is the current image rendition choice.

Icons	Image	T O
R -	Rescale Image	
lp	12 Level Gray	
	▶ 16 Level Gray	
MC	12 Color	
9		
499	▶ Show Baseline	
01	Hide Baseline	

While a menu is visible, use the mouse to move the pointer. To choose a command, highlight the line you want and click. To clear the menu from the screen, move the pointer into a clear area of the display and press the mouse button. The last menu you displayed disappears.

Certain XR-26 users may find it distracting when menu's drop down whenever the pointer is moved across the top of the display. XR-26 software permits the operator to choose between "Click" and "No Click" drop-down menu operation. Please refer to Section 5.4.5 for further information.

1.4.3 KEYBOARD EQUIVALENTS TO MOUSE TECHNIQUES

<u>Technique</u>	<u>Keyboard Key</u>
Click	Press and release the HOME key once.
Double-click	Press and release the HOME key twice, quickly.
Drag	<ol style="list-style-type: none">1. Press and release END key.2. Use the arrow keys to move pointer where desired.3. When finished dragging, press HOME key.
Shift-click	<ol style="list-style-type: none">1. Use arrow keys to place pointer on the first object you want to select.2. Press and release the HOME key to select the object.3. Use arrow keys to move pointer on top of next symbol you want to select.4. Press SHIFT and HOME keys to select new object.5. Repeat steps 3 and 4 to select any additional objects.

Table III - Keyboard Equivalents to Mouse Techniques

In the instance where an XR-26 operator does not have a mouse available to manipulate the pointer, it is a simple matter to use appropriate keys on the computer keyboard to simulate the mouse commands. Table III lists each mouse command next to its equivalent keyboard command.

NOTE: If the arrow keys don't move the pointer, press the CTRL key and try again. The CTRL key switches the keyboard into and out of the pointer movement mode.

As an operating convenience, XR-26 software will automatically perform the command whose button is displayed with a heavy dark border whenever the **Enter** key is pressed. This heavy-bordered button is referred to as the Default button. For example, many setup screens display the **Accept Changes** button. Notice that this button has a thicker border than the other button. Pressing the keyboard's **Enter** key causes the **Accept Changes** command to be performed, then the next logical screen appears. This convention is true regardless of the presence of a mouse. Should you be editing some on-screen dialog and you press the **Enter** key, the command shown in the thick-bordered button is performed and the next logical screen appears. If you weren't finished editing dialog, you need to return to the previous screen.

The numeric keypad segment shown in Figure 2 indicates the HOME and END keys, as well as the arrow keys. Make sure the "Num Lock" indicator on the keyboard is not illuminated, or the commands listed in Table III will not work.

Some keyboards provide duplicate HOME, END, CTRL, and other control and cursor movement keys. When reviewing the keyboard-equivalent commands, be aware that these duplicate keys perform functions identical to numeric keypad keys.

7 Home	8 ↑	9 PgUp
4 ←	5	6 →
1 End	2 ↓	3 PgDn

Figure 2 - Numeric Keypad

1.4.4 DIALOG EDITING TECHNIQUES

Much of the XR-26 program is controlled with the mouse and drop-down menu interface, making only limited use of the computer keyboard. However, the keyboard remains the principal method of entering patient personal data, scan comments, and other information. This data is entered using Dialog Editing Techniques.

Pressing the numeric keypad or other keyboard key named below at the left causes the action described to its right. The [§] symbol indicates that the effect is continuous if the key is held down.

<u>Keyboard Key</u>	<u>Edit Cursor Action</u>
4 - ←	Moves the text cursor one character left without deleting. [§]
6 - →	Moves the text cursor one character right without deleting. [§]
2 - ↓	Moves the text cursor to the next field in which you can enter or change information. [§]
→ (Tab)	Moves the text cursor to the next field in which you can enter or change information. [§]
8 - ↑	Moves the text cursor to the previous field in which you can enter or change information. [§]
Del	Deletes character to the right of the text cursor. [§]
←	Deletes character to the left of the text cursor. [§]
Esc	Erases all characters in the field where the text cursor is located.

To select any field shown to enter data, use the mouse to position the pointer in that field, then click. A narrow vertical line appears just to the right of any information present in that field. Use any alphabetic or numeric symbol on the keyboard to enter the information you desire. Then either click on a new field for data entry or press the TAB key to move to the next available field. If you hold down the SHIFT key as you tap the TAB key, the narrow vertical line is moved to the previous field.


Remember that if the **Enter** key on the keyboard is pressed while in a dialog edit screen, the command appearing in the thick-bordered button is performed, then the next logical screen appears. If you weren't finished editing, you need to return to the screen you were working on to complete your input.

Dialog Edit Example

The following is an example dialog edit session used for entering a scan comment. The techniques described may be used for any dialog entry you care to do.



First, select **Analysis** and click on the **Edit Scan Comments** command.

On the Edit Scan Comments screen shown at the right, notice a thin line at far left of the top comment line. This symbol is the character entry cursor. It is automatically positioned at the end of newly-entered text and moves as the text is received.

Use the computer keyboard alphabet and numeral keys to enter comment text as desired on the first line. To move to the next comment line, either press the **Tab** key or use the mouse to click on the next line.

To make corrections, click on the inaccurate line, then use the appropriate arrow movement key to move the line symbol next to the error. Make the correction(s) as needed.

The completed comment screen is seen at the left. This information will appear on any printed reports of this scan. To accept the final text entry, click on the **Accept Changes** button.

This page reserved.

2. DAILY OPERATIONS

2.1 ROUTINE OPERATIONS

Before scanning a patient, it is necessary to turn-on the various components of the instrument, check that each is operating properly, and to have on-hand some pre-formatted data diskettes. Note that these diskettes should NOT be write-protected. Please refer to Section 5.3.5 for formatting instruction.

2.1.1 POWERING UP THE SYSTEM

Check the Laser Manual Shutter paddle (found below the scanner arm - see Figure 4, later in this section) to make sure that it is in its horizontal position. Lifting the paddle to its horizontal position opens the shutter and permits the laser to project onto the scanner table.

Turn on the power to the scanner by depressing the rocker switch located in the power cord receptacle at the lower right side of the instrument. The switch is rather small and is positioned at the left-most side of the power receptacle. Norland does not recommend any warm-up period before calibration or patient scanning.

Next, power-up the Control/Analysis work station: Turn on the computer display by pressing the switch to the ON (the 'I') position. Next, turn on the computer by flipping the power switch to the ON (the 'I') position.

A series of characters will appear in the upper left corner of the computer screen. This is an indication that the computer is doing its self test. Soon, the 'C>' prompt appears. The 'C>' prompt may actually resemble 'C:\>' or something else. Regardless of appearance, 'C>' will be used herein to denote the root DOS prompt.

To begin scanning and analysis operations, type 'XR-26' at the 'C>' prompt, then press **Enter**. A screen similar to the Standard Idle Display for XR-26 Version 2 Software appears, as shown in Figure 3.


Figure 3 - Main Menu or Standard Idle Display screen

In the event of a start-up problem or other difficulty, please refer to Section 5.9.4 for assistance. If further questions arise, please refer to training material that may have been given to you at XR-26 installation, or discuss them with the Norland Customer Service technician.

Idle / Main Menu Display

The Idle display (Figure 3) includes Norland's address and telephone number and the revision level of the software installed in the computer and the scanner. The XR-26 system automatically completes initialization procedures. If the word "None" appears next to the 'Scanner Software Revision:' line, check to be sure the cable between the computer and scanner is properly connected. Refer to Section 5.9.4 if necessary.

Complete the power-up sequence by turning the printer on. You are now ready to make selections from the Main Menu bar titles to begin setup or analysis operations. The instrument does not require any warm-up period before scanning operations.

XR-26 Software

As explained earlier, Version 2.0 and later XR-26 software is devised within the pictorial, visually-oriented structure of the GEM®/3 operating environment. XR-26 software is an image-based program designed to use to full advantage the excellent graphics and information-handling capabilities of GEM®/3. XR-26 software is also a 'mouse'-based application. All program actions are initiated by moving an arrow-shaped pointer around on the screen, as described in Section 1.4.1. Please review those instructions since the operator interface terminology defined there is used exclusively throughout this guide.

One of the most important aspects of using a personal computer is learning how to preserve accumulated data. Although safeguards have been built into the software to preserve the integrity of acquired Patient Data, it is possible that a mechanical failure could destroy data. For this reason, Norland strongly recommends that a routine hard disk backup procedure be established. On a daily or weekly basis, the critical XR-26 patient file data should be backed-up on diskettes, thus reducing the chances of significant data loss. Please refer to Section 5.3.5 for further information.

In review, there are a total of five items that need to be turned on, as follows:

1. Scanner Rocker Switch
2. Laser Manual Shutter Paddle
3. Computer Monitor
4. Computer
5. Printer

2.1.2 POWERING DOWN THE SYSTEM

The scanner should be turned off when not in use for periods longer than two hours. This will lengthen the service life of the x-ray source components. There is no warm-up period required before performing calibration or patient scans.

The order in which the equipment is turned OFF is not critical. To leave the XR-26 program, select **File** and click on the **Quit XR-26** command. Click on the **Exit to DOS** button to exit to the 'C>' prompt. Another way to exit XR-26 program is to click on the **Close Box** symbol (described in Section 1.4.), then click on the **Exit to DOS** command.

If printouts have not been completed, a message appears asking if you wish to cancel the remaining printouts. If you want to complete the printouts, wait until they are finished before exiting the XR-26 program. To cancel printouts and discard the print files, click on the **Terminate** command, then click on the **Exit to DOS** command. The printouts may be requested again the next time you use XR-26.

At the 'C>' prompt, you may turn off the computer. If your facility has a policy of leaving computers on at all times, it is acceptable to leave the computer on.

2.2 SCANNER CONTROLS AND INDICATORS

All of the scanner controls and indicators are described in this section. Figure 4 depicts the various scanner control elements.


Figure 4 - Location of Scanner Controls

CAUTION: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

The Laser Manual Shutoff knob is located at the lower distal end of the scanner arm on certain XR-26 instruments. The Laser Manual Shutoff knob may be rotated fully counterclockwise to avoid endangering the

patient's eyes. The Laser Manual Shutoff paddle is located on the under-side of the scanner arm on other XR-26 instruments. If the patient's eyes are jeopardized by the laser beam, either slap the paddle to it's closed (down) position or press the LASER OFF button on the Scanner Control Panel.

Two control areas exist on the scanner. The Scanner Control panel is located at the distal end of the scanner arm, and is depicted in Figure 5. The other control area is at the lower right end of the scanner pedestal. The combination power switch/fuse receptacle is located here, as is the communication cable connector. The scanner patient pad is easily lifted off the scanner to access the x-ray source and other inner mechanisms. For identification purposes, the power switch end of the scanner pedestal is designated the 'head' of the table. The other end is called the 'foot' and is the location of the 'Table Origin' or Home position.


Figure 5 - Scanner Control Panel

The Scanner Control panel includes several touch-sensitive keys and operation indicators. The **MARK** key is pressed when the operator is asked to mark an extent of the scan site. The **EMERGENCY OFF**

key is used to instantly stop all scanner operations and turn off the x-ray beam and laser. Separate indicators are illuminated whenever the X-RAY source has power applied to it, whenever the SHUTTER is open, and whenever the LASER has power applied to it.

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength. **DO NOT STARE INTO BEAM.**

The laser is turned ON by pressing the control panel key labeled **LASER ON/OFF**. The laser is only enabled when its operation is needed for the procedure being performed. After turning the laser on, a red dot appears directly below the scanner arm. If the red dot is not visible or not perfectly round, check the Laser Manual Shutter paddle to make sure that it is fully elevated (horizontal).

Please note that if the **EMERGENCY OFF** key is pressed, the scanner stops all movement, closes the shutter, and turns off the laser and x-ray power. The scanner's power switch must be turned off and then on again to recover from this condition. Scan procedures may then be re-initiated once the reason for invoking the **EMERGENCY OFF** condition has been resolved.

Also included on the scanner control panel are the scanner arm movement keys. The arrows indicate the direction the arm and source will move if that key is pressed. Arm movement commands are requested by screen messages when appropriate.

2.3 QUALITY CONTROL

Quality Control consists of the procedures used to assess and maintain the XR-26's measurement competence. Faithful performance of the techniques outlined in this section will enhance the operator's level of confidence in the actual bone mineral density measurements. Quality Control includes calibration, management reports, and phantom scanning.

2.3.1 CALIBRATION

Calibration is the process of adjusting the XR-26 instrument to compensate for spectral shift in the x-ray source, slight differences in the interface circuitry, and other measurement factors.

Two XR-26 calibration processes have been implemented. An initial 12 hour 'characterization' scan is needed to accumulate enough statistics from the scanner to fully initialize measurement values. This procedure is usually performed at the time the instrument is first installed at the site. However, replacement of a major component of the XR-26 scanner may dictate another characterization scan. Instructions for performing this procedure will be supplied by Norland in the event such re-characterization is required.

The second calibration process should be performed daily and is completely automatic. This daily calibration task takes about 20 minutes, from initiation through completion of printout. To begin calibration, select the **Scanner** command and click on the **Begin Calibration** command.

If the **Begin Calibration** command is de-emphasized (not able to be selected), you will need to enter the Calibration Standard Identification

values as described in Section 5.4.1. Ordinarily, this information is entered by the Norland Technician at the time the instrument is installed. Once entered, the *Begin Calibration* command is activated and is used to initiate the calibration process.

The operator is asked to place the calibration standard in its position on the scanner table. The bracket on the calibration standard incorporates keyhole-shaped openings. Simply slide the openings over the heads of the mounting posts and permit the standard to drop to its resting position. The mounting posts are at the left (foot) end of the table. The scanner arm should be positioned near the foot of the table.

Click on the **Cancel** command if the procedure needs to be terminated for any reason. Please note that if the **EMERGENCY OFF** key is pressed, the scanner's power switch must be turned off, then on again, to reset the unit. As the screen message requests, switch **ON** the laser beam.

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength. **DO NOT STARE INTO BEAM.**

Move the laser spot as close as possible to the 'A' point inscribed on the top of the standard, using the scanner arm movement keys. When satisfied with the laser dot's position, press the **MARK** key. Should you wish to cancel the procedure, click on the **Cancel** command.

After point 'A' is marked, a message asking that the 'B' point on the calibration standard be **MARKed**. Locate the 'B' point with the laser beam and press the **MARK** key. To re-mark point A, click on the **Remark Point A** button.

When point 'B' is marked, the XR-26 turns the laser off and completes the calibration procedure without further operator assistance.

A series of instrument diagnostics are next performed, indicating pass or fail status for reference beam, shutter, high, medium, and low

attenuation filter tests. If an error is detected, a message explaining that calibration is not able to be completed appears on the screen. It is advised that a written note be made of the error codes, then Norland Customer Service should be notified.

After the quality control diagnostics have been successfully completed, a scan of the calibration standard begins. As the scanner arm moves to its pre-programmed positions and measurements are made, the display screen is updated. The estimated time remaining for the procedure is displayed. Calibration results appear when the process is complete. The results of the successful calibration are saved in disk memory.

Clicking on the **Cancel** button any time before the Calibration Results screen appears causes the procedure to be terminated, and the calibration and quality control measurements made up to that point are discarded. The instrument will then use the previously-recorded values.

At the end of the Calibration procedure, a complete summary of calibration measurements is presented. To obtain a hard copy of the calibration results, select **File** and click on the **Print Report** command.

2.3.2 PHANTOM SCAN

Each XR-26 instrument is supplied with a Norland Anthropomorphic Spine Phantom. This accessory is constructed of calcium sulfate and cast polyester. The phantom is fully characterized at the factory and the resultant values are recorded on a label attached to it.

As a quality control index, the phantom should be measured once after each daily XR-26 calibration scan using a 12 cm Scan Width and techniques described in Section 3.1, Spine Measurement. Single (daily) measurements of L2-L3-L4 BMD, BMC, and Area, should fall within 2.0% of the factory assigned values.

Routine daily scanning of the phantom provides a 'rough' accuracy validation, as well as a long-term precision check. To more closely

monitor short term system performance, occasionally measure the phantom several times in a row (not less than six). The L2-L3-L4 mean values for BMD, BMC, and Area should fall within 1.0% of the factory characterization values.

The coefficient of variation (CV) of these results should be less than or equal to 1.0 percent, where:

$$CV = \frac{(\text{standard deviation} \times 100)}{\text{mean}}$$

Please refer to Section 5,6 for further detail about this CV specification.

Verifying the Area measurement serves to verify the XR-26 step size (mm per line) because one of the Area parameters is the number of scans lines required to cover the length of the phantom. The XR-26 software converts the number of scan lines into centimeters (cm) when performing the Area calculation (cm²). This verification is important because serial studies dictate the use of the same anatomical lengths from scan image to scan image.

If measurement values exceed the specified range, the measurement technique should be reviewed and the phantom re-scanned. If the out-of-specification condition persists, Norland Customer Service should be notified.

This completes the descriptions of the Scanner Controls and Quality Assurance subjects. We invite you to review other sections of this Operator's Guide.

This page reserved.

3. A-P SPINE MEASUREMENT

This spine measurement mode in the XR-26 is designed specifically for bone mineral quantification of the lumbar spine, using an anterior-posterior projection. The region of interest is the L2-L3-L4 segment, which is analyzed both in total and as individual vertebrae. The analysis software excludes vertebral transverse process areas from the bone mineral calculations. It also provides a means of excluding scan artifacts. A Lateral Spine Measurement Package is available as an option for the XR-26.

Proper operation of the XR-26 and interpretation of its measurements requires a thorough understanding of the factors affecting precision and accuracy. Such information is provided in Section 5.6, "Precision and Accuracy."

Although a female patient is used for procedure descriptions herein, the same techniques are to be used with a male patient. In addition, references are made to various sections in the Base Operator's Guide. You should become familiar with this referred material before performing the techniques described herein.

The scanner is equipped with a low power helium-neon laser which is used as a positioning aid. The red laser light spot is projected onto the patient and indicates where the x-ray beam will pass.

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength.
DO NOT STARE INTO BEAM.

3.1 SCANNING THE SPINE

3.1.1 Preparing Patient

The patient may wear a cloth garment. There should not be any metal or high-density plastic objects in the scanning area that could interfere with the scan, especially back braces, long chains or other jewelry; any metal on their clothing such as rivets, buttons, zippers, or other types of fasteners. In other words, check for any item that would add a density to the study. If any of these items are found, be certain to keep them clear of the area to be scanned.

The leg rest, a large fabric-covered block, should be readily available. The different dimensions of the leg rest accommodate various patient leg lengths.


Figure 6 - Patient on Scanner with Leg Block

The patient lies down in the center of the patient pad with her right side adjacent to the scanner backrest (this places her head at the 'head' end of the scanner). The scanner arm should be positioned over the patient's midsection using the scanner arm movement keys. The patient is then arranged so that her legs are comfortably positioned on the leg rest, as shown in Figure 6. The use of the leg rest reduces lordosis of

the spine and helps the patient remain comfortable. Also pay attention to patient placement so that the scanner arm doesn't touch the patient. It is good practice to develop a consistent positioning method for all patients.

As the patient may be lying in this position for several minutes, it is important that she is made as comfortable as possible. This care reduces the chances of extraneous movement. Make sure she is warm or cool enough, and that there is a pillow for her head. The use of a sheet or light blanket for warmth or modesty will not interfere with the scan results. The patient is now ready for scanning.

Select Patient

Please refer to Section 5.1.1 for detailed instruction on using the patient selection features of the XR-26. If the patient to be scanned is a new patient and needs to be added to the patient directory, please refer to 5.1.1.

3.1.2 Scan Preparation

After selecting the patient, select **Scan** and click on the **AP Spine** command. A screen requesting the patient's height and weight, and the names of the scanning technician and physician appears. This information is entered using Dialog Edit Techniques described in 1.4.4 and is automatically added to the patient's data file. This information is requested each time a patient is scanned.

With the patient in position, the operator should prepare the scanner as prompted by the display. Turn on the laser and observe that a red dot appears directly below the scanner arm. If the red dot is not visible, check that the Laser Manual Shutter paddle is fully elevated (horizontal).

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength. **DO NOT STARE INTO BEAM.**


Figure 7 - Recommended Scan Area

The typical lumbar scan region is shown in Figure 7. Using the umbilicus and xiphoid process as center line references, the scan extends 7 cm to either side of center and from mid-L1 to just beyond the iliac crests. To include all of L2 through L4 in the scan, the start point should be marked about 3 cm below the xiphoid process. When just beginning doing scans, however, it is usually helpful to start at the xiphoid process so that the twelfth rib appears in the scan.

Position the laser dot over the scan start point and press the **MARK** key. A beep is heard, indicating that the computer acknowledges the position. The screen message now asks that the end point be identified. If you click on the **Cancel Scan** button, another message appears, asking if you wish to terminate the operation. Click on the

Resume button to continue or on the **Terminate** button if you wish to cancel the scan.

It is helpful to include part of the iliac crest in the scan image. Including this feature makes identification of L2, L3, and L4 much easier. Move the laser dot down through the umbilicus on the patient's center line, about 3cm below where the beginning of the iliac crest can be felt. Take care to avoid bumping the patient's raised legs with the scanner arm. Press the **MARK** key. The ensuing beep indicates that the computer acknowledged the scan end point. The initial and final extents so **MARK**ed have defined the length of the scan area.

If you click on the **Remark Start Point** button, the computer re-displays the Initial Mark information box. Re-position the scanner arm at the desired start point and press the **MARK** key to re-mark the initial point. Otherwise, click on the **Cancel Scan** button to terminate the process.

Once the initial and final points have been marked, a Scan Parameters screen similar to Figure 8 appears, listing point resolution, line resolution, scan width, and scan speed parameters. Point and line resolution and Scan Speed parameters are fixed by XR-26 software and cannot be changed. 13.95 cm is the Scan Width recommended to provide ample soft tissue characterization. If desired, the scan width may be redefined using the Dialog Edit Techniques, Section 1.4.4. However, subsequent scans of the same patient should use the same Scan Width.


Figure 8 - Scan Parameters screen

If necessary, click on the **Remark Scan** button. Instructions for scanner arm placement appear. Follow the on-screen guidance to re-mark the scan. Or, if you click on the **Cancel Scan** button, you are asked to confirm your desire to terminate the operation. Click on the **Resume** button to continue or on the **Terminate** button if you wish to cancel the scan.

3.1.3 The Spine Scan

After viewing the scan parameters, click on the **CONTINUE** button. An X-Ray Parameters screen similar to Figure 9 appears, indicating the intended Anode Voltage (in kilovolts), Anode Current (in milliamps), and Estimated Scan Time (in minutes).


Figure 9 - X-Ray Parameters Information Box

If you concur with the scan time, caution the patient to lie still, then click on the **Start Scan** button. The XR-26 turns-off the laser and selects the filter combination dictated by patient thickness, then measures a baseline. The computer may display a message stating that the x-ray source is being brought to operating condition. The scanner arm will soon begin to move back and forth across the patient. The system completes the entire scan now without intervention. An audible beep is emitted by the computer at the end of a complete scan.

If you do not concur with the listed scan time, click on the **Edit Parameters** button. Make the changes you desire to Scan Width or Scan Speed, and click on the **Continue** button. The computer will display the effect of this change on the Scan Time shown in the ensuing X-Ray Parameters screen.

The Scan Progress screen shows how many lines will be scanned and gives an estimate of scan time remaining. While the scanner is acquiring data, an image of the scan is built up on the screen. It is important for the operator to observe the screen image as it appears because the patient may shift position. Close observation permits the scan to be terminated, if necessary, thereby minimizing patient radiation exposure.

The image is rendered in a gray scale (or color spectra, if selected). Note the Density Spectrum that appears near the image. The 'L' indicates the lowest detected bone density and the 'H' indicates the highest bone density. Intermediate bone densities are evenly allocated to the intermediate gray (or color) levels. These images aid in placing

the cursors and are not intended for diagnostic purposes. Subsequent analysis is performed on the actual data acquired.

Remember that if the **EMERGENCY OFF** key on the Scanner Control Panel is pressed, the scanner must be turned off, then on again. Scanning may then be resumed.

Once scanning has started and the **SCAN PROGRESS** screen is displayed, clicking on the **Cancel Scan** button causes the current scan line to be completed. An information box appears, asking if you wish to terminate the operation. Click on the **Resume** button to continue.

If you click on the **Terminate** button, a **Scan Terminated** screen appears (very similar to Figure 10). You are asked if you wish to **Analyze**, **Rescan**, or **Discard** the data. Examine the displayed image and determine if it suits your purposes. If not, click on the **Rescan** button to perform a rescan, as described in Section 3.1.4. To discard the data, click on the **Discard** button and confirm the discard operation. Clicking on the **Analyze** button permits analysis of the acquired data, as described in Section 3.2.


Figure 10 - End of A-P Spine Scan screen

When the scan is completed normally, the XR-26 presents the **Scan Complete** screen, as seen in Figure 10. You are asked if you wish to **Analyze**, **Rescan**, or **Discard** the data. Click on the **Analyze** button to continue the procedure, as described in Section 3.2. If you click on the **Rescan** button, you may perform a rescan of the patient, as described in Section 3.1.4. If you click on the **Discard** button, the data is not saved and the **Main Menu / Idle** screen eventually reappears.

3.1.4 Rescan Patient

When the scan is completed, the operator normally clicks on the **Analyze** button, then follows the procedures described in Section 3.2. However, if the scan image indicates that a rescan is needed, the operator may click on the **Rescan** button. A reduced-size scan image appears with two cursors superimposed. This image provides an idea of the underlying anatomy or the existence of scan artifacts.


Figure 11 - Rescan setup screen

The "X" cursors are to be positioned on the image to define the start and end mark points for the rescan. Drag the upper or lower "X" cursor to where you want it on the screen image. Then drag the other "X" cursor to the desired position.

Notice that as you move the "X" cursors, a rectangular region showing the 'scan window' moves with it. This feature, seen in Figure 11, should greatly assist placement of cursors to exactly define the scan region. The "X" cursors are connected by a line that bisects the 'scan window' and indicates the axis of arm movement. It should be used to guide cursor placement so that the ensuing scan is centered on the patient's spine.

Caution the patient to remain lying quietly, then click on the **Edit Scan Parameters** button. Verify the scan width value shown on the Scan Parameters screen. Click on the **Continue** button, then on the **Start Scan** button. The scanner arm will begin to move across the patient and the scan will be completed without further intervention.

3.2 SPINE ANALYSIS

The information presented in this section continues with the normal scan-then-analyze operation sequence. XR-26 software provides a way to re-analyze scan data, as described at paragraph 3.2.4, below.

In a Spine Analysis, the characteristics of three lumbar vertebrae (L2, L3, and L4) are calculated based on the information stored in the scan data file. Data analysis uses the bone and soft tissue values from the scan data to produce an image of the scanned anatomy on the display. To proceed to Analysis, click on the **Analyze** button.

The analysis process begins with an initial vertebral definition. The software makes additional refinements to the calculations, then the cursor positions are reviewed and verified. If scan artifacts need to be removed from analysis, an Include/Exclude process is provided. Finally, numeric results are calculated and presented.

If you click on the Main Menu bar **Analyze** title after making analysis changes, you are asked to *Save*, *Abandon*, or *Cancel*. *Save* means save the data to disk and return to the Idle screen. *Abandon* means to abandon the work done so far, discard the data, and return to the Idle screen. *Cancel* means to return to the screen at which the Main Menu command was clicked on. If you are following this explanation, you should click on the **Cancel** button to continue.

3.2.1 Initial Vertebral Definition

The Initial Cursor Positioning display (Figure 12) shows the various tissue densities found during the scan in a color or grayscale density spectrum. Please note that the cursor positions illustrated here are as

they originally appear and do not represent recommended cursor placement. They must be repositioned as described below. Also note that not enough of L2 was acquired, so a rescan would be recommended if this were an actual patient scan.


Figure 12 - Initial Cursor Positioning screen

The analysis process begins with a preliminary placement of the intervertebral cursors. Four line cursors extend across the width of the screen and have been placed in pre-set positions. These cursors are to be positioned by the operator so that they cut the lumbar spine image at appropriate intervertebral spaces, thus defining the L2, L3, and L4 vertebrae.

To the right of the spine image is a bone mineral content histogram. The BMC histogram is a graphical presentation on a scan line-by-scan line basis of the vertebral bone mineral content. Histogram lines which are short denote less bone mineral content, longer histogram lines denote higher bone mineral content. In most scans, the histogram, plus the actual image, will be all that is needed to identify the vertebrae. In the normal spine, the BMC value will be lowest at the intervertebral gaps, so simply set the cursors at these minima in the BMC histogram.

When analyzing subsequent spine scans of the same patient, it is important to be able to position the intervertebral cursor group as a unit without changing their size or positional relationship since the cursor placement and size must be the same as that used for the initial scan. This is accomplished by a shift-click anywhere on the image itself (not on the histogram). The description 'L2-L4' appears under the "Current Region" legend, but the width and height of this cursor is not specified.

Each of the four cursors is adjusted for appropriate intervertebral separation. Click on the small box centered on the bottommost cursor.

This box will be referred to as a Control Point, Observe that "L4/L5" appears beneath the Current Region legend. Drag it into position at the bottom of L4. Notice that as the cursor is moved, the value displayed for L4's length changes. Click on the **Next Region** button or on the Control Point for the next cursor up and drag it into position between the bottom of L3 and the top of L4. Similarly, position the next cursor between L2 and L3, and the last cursor at the top of L2.

When satisfied with cursor placement, click on the **Continue** button. The computer performs additional calculations, then presents the Review Regions screen (Figure 13).

3.2.2 Review Vertebral Definition


Figure 13 - Review Regions (cursor placement) screen

Verify that the positions of the cursors (as depicted in Figure 13) are appropriate. The cursor positions so defined are stored with the analysis.

Baseline data (non-bone points) are shown or 'hidden' by clicking on the appropriate **Image** command. Removal of baseline data reveals all scanned points corresponding to the chosen display spectrum

without regard to their 'bone' or 'non-bone' status. This feature may be helpful in locating anatomical landmarks such as ribs and iliac crests, which have been automatically excluded from the analysis. To select an alternate image rendition, select **Image** and click on the rendition command of your choice.

If artifacts have been included as a part of the analysis, they may be removed using the INCLUDE / EXCLUDE feature, as described in Section 3.2.5.

At any subsequent spine scan and analysis of the same patient, the XR-26 software will automatically recall the cursor placement, maintaining the same vertebral lengths. These cursors should be positioned as a group to define the same analysis regions as were used for the initial scan/analysis. In particular, the vertebral lengths must be the same.

Direct comparison with a printout of a previous scan helps in accurate repositioning of the intervertebral cursors. Another method of positioning for a second (serial) scan of a patient is to use the Comparison Image feature. The Comparison Image is derived from the analysis of the first scan performed on the subject patient and is presented in the center of the Define Regions display.

Figure 14 shows an example side-by-side Comparison Image. To obtain a similar Comparison image, select *Image* and click on the *Show Comparison* command. The analyzed image from the first scan of this type for the patient is recalled and presented in the center of the display, using the same linear scale as the existing image. In the case of the AP Spine scan, the histogram will be replaced with the Comparison Image.


Figure 14 - Show Comparison image with original scan at the right.

In summary, the XR-26 aids the operator by placing the cursors using the same spacing and positioning as the initial scan; by moving the cursors together as a group; by displaying the segment length (cm) on the screen (and on printouts) for reference; and by providing a

Comparison Image. Improper cursor positioning may render the analysis results inappropriate for serial studies.

NOTE: When performing follow-up scans on patients, it is very important to be consistent; that is, to re-position the intervertebral cursors as close as possible to the original vertebral placement. It is also important to use the same Scan Width as was originally used (13.95 cm is recommended).

3.2.3 Numeric Results

When cursor positions have been verified, select **Results** and click on the **Summary Results** command. The computer performs summary calculations and stores the results. An example Summary Results display is shown in Figure 15.


Figure 15 - A-P Spine Scan Summary Results screen

The image is not intended for diagnostic purposes. The image rendition or baseline selection may be changed by selecting **Image** and clicking on the appropriate command.

On the Summary Results screen, the intervertebral cursors appear in position over the image. Listed are measurement results for each region and the operator-selected vertebral lengths (for reference), such as bone mineral density (BMD) and bone mineral content (BMC). An estimate of the coefficient of variation (CV) for each measurement is included in parentheses beneath it.


Figure 16 - Detailed Results Screen, Page 1

	BMD	BMC	AREA	LENGTH
	g/cm ³	g	cm ²	cm
L2	0.823	8.72	10.59	3.15
L3	0.802	11.47	14.31	3.45
L4	0.799	11.23	14.23	3.15
L2-L4	0.802	31.43	39.13	9.75

Figure 17 - Detailed Results Screen, Page 2

XR-26 Software provides two additional Results screens. The Detailed Results page 1 screen, seen in Figure 16, presents the scan image alongside the patient personal data. The Detailed Results page 2 screen lists the values measured for each vertebral region, area, and total vertebral column measurement values (see Figure 17).

Comments about this scan, such as special conditions or observations, may be entered using Dialog Edit Techniques (1.4.4) after selecting **File** and clicking on the **Edit Comments** command. Any comment so entered will be saved with the data and is included on any subsequent numeric results printout. Click on the **Accept Changes** button to preserve the scan comments.

EDIT SCAN COMMENTS

This is the Scan Comment. The _____ technician or doctor may note _____ exceptional or unusual conditions _____ for future reference. _____

3.2.4 Reviewing Previous Results

When an analysis for a particular patient is desired, it is necessary to specify the patient and the scan of interest. To do so, select **File** and click on the **Select a Patient** command. Double-click on the particular Patient name, then double-click on the desired spine scan data file. The computer then retrieves the data for the selected scan and

presents a Results screen. More detail about the patient or scan selection features of XR-26 software may be found at Section 5.1.1.

If the data was stored on diskette, a "Loading Data" message requesting the insertion of that specific Volume disk is displayed. Place the diskette in the disk drive and click on the **Continue** button.

For Summary and page 1 Detailed Results, intervertebral cursor positions are recalled with the data and are displayed in position. Page 2 Detailed Results lists measured values for each region.

If cursor positioning needs to be redefined, select **Regions** and click on the **Modify Regions** command. The Review Regions display is presented and the desired changes may be made, as described above (paragraph 3.2.2). To view a Comparison Image, select **Image** and click on the **Show Comparison** command. If additional artifacts have been identified, select **Regions** and click on the **Include/Exclude** command. Mark and process the artifacts as described below.

3.2.5 Artifact Include / Exclude

The spine analysis program has been designed to automatically exclude any 'bone points' which are not connected to the spine. Similarly, it will automatically eliminate any voids or islands of soft tissue surrounded by vertebral bone. Occasionally, however, it may be necessary to manually remove non-spine artifacts, such as ribs, iliac crests, calcified aortas or extra-osseous calcifications, previous callous formation, etc., from the measurement region by using the Exclude function.

To exclude such artifacts, from the displayed Results screen, select **Regions** and click on the **Modify Regions** command. Reselect **Regions** and click on the **Include/Exclude** command.

On the Include/Exclude screen, a small box appears over the image (as in Figure 18). It will be used to circumscribe such artifacts for removal

from analysis consideration. The entire cursor box may be moved up, down, left, or right.

Position the pointer anywhere inside the box and drag the box to the position desired. Then, drag an edge of the box where desired. Continue with this positioning sequence for the remaining sides of the cursor.

When satisfied with the box size and placement, click on the **Exclude** button to exclude the region from results consideration.


Figure 18 - Include/Exclude screen

If you wish to include a portion or all of an excluded region, especially islands of tissue, click on and drag each edge of the cursor box to circumscribe the area of interest. When satisfied with cursor size and placement, click on the **Include** button. Data from the newly-defined region is included in any subsequent results report.

This process of excluding each fragment of the artifact may be repeated as often as needed. When satisfied with region exclusion, click on the **Continue** button. Review the intervertebral cursor positions, then select **Results** and click on the desired results command. Data from any excluded region(s) is removed from the tabulated results. Data from any included region(s) is added to the tabulated results.

3.2.6 Alternate Analysis Technique

In some patients, the spine has deteriorated to such an extent that it is impossible to locate one or more intervertebral spaces. Techniques have been developed that may be used to follow the patient serially. One of these techniques is to use a fixed-length region.

In order for the fixed-length region to be accurately re-positioned, the scan images must include both iliac crests. To set up a fixed-length region, the single-line cursor that is normally positioned at the bottom of L4 should be positioned across the top of the iliac crests. It may be best to hide the baseline so all scanned points may be seen. The single-line cursor that is normally positioned at the top of L2 should be positioned some fixed-length above the iliac crests. For example, positioning the L2 (top) cursor 70 lines above the iliac crests results in a fixed length segment of 10.5cm (70 lines x 1.5mm per line). The middle two single-line cursors should be either ignored or moved to convenient locations between the top and bottom cursors. Once the cursors are positioned, the operator may continue with the analysis. The values given in the Results display for these two cursors should be ignored. Only the Total values for the L2-L4 segment should be used for comparison purposes.

If some of the vertebral separations are able to be determined, place the appropriate single-line cursors at those separations. The remaining single-line cursors (those not able to be positioned) should be placed anywhere convenient. Only the values for the properly-positioned cursors should be used. It may be useful to note this cursor placement in the Scan Comments record.

3.2.7 Rescale Image and Reanalyze Scan

It is possible that a scanned region includes a very high density object. The Report image density spectrum is scaled according to the range of densities found in the scan. An high-density object included in the analysis may cause the resultant density spectra to be spread across a wider-than-normal range of density values, so the on-screen and printed density rendition may be skewed. The portrayal of lower-density scan elements may be obscured because of this skewed density spectrum.

XR-26 software permits these high-density objects to be removed from analysis consideration using the Exclude feature. Once removed, the scan image may be rescaled so that a more normal range of density values is portrayed. To do so, use the instructions given above to

exclude the object. Then select *Image* and click on the *Rescale Image* command. The software will disregard the excluded areas from the analysis and present the image in a more typical density spectrum.

The Include feature may be used to recover each excluded portion of a scan in a piece-by-piece fashion. The Reanalyze operation may be used to recover all excluded regions in one operation, without moving any intervertebral or Special Region cursor(s). To reanalyze a scan and recover all excluded regions, select *Regions* and click on the *Reanalyze* command. Confirmation of this operation is NOT made, so be sure you have no need for any of the existing excluded regions. The analysis process starts over from the beginning.

The High Density Detection feature can discard any data point containing a very high density value, thereby removing these artifacts from the computed measurement values. It is useful in situations where an implanted prosthetic device is in the scan region, such as bone screws and fixtures. High Density Detection may be enabled or disabled at the Analysis Results screen by selecting *Regions* and clicking on the desired *High Density Detect* command. Once set, this mode remains in effect until changed.

3.2.8 Printed Report

The Summary Report, shown in Figure 19, depicts essentially the same information seen at the Summary Results screen. It takes less time to generate and print than the Detailed Report, which combines the information seen in both Detailed Results screens on one piece of paper. The printed image is not intended for diagnostic purposes.

Choices made in the *Image* Menu are used when generating hard-copy printout of the analysis. These choices include images rendered with or without baseline data and in color or gray scale. What is seen on the display is what will be printed. The number of copies that will be printed is pre-selected at Printer Setup, described in Section 5.4.4.


Figure 19 - A-P Spine Analysis Printed Report

If a Summary report is desired, from the Summary Results screen, select **File** and click on the **Print Report** command. Then click on either the **Print Immediately** button or **Add to Queue** button. The preselected number and rendition of the analysis results is obtained. If additional printed reports are desired, click on the **Print Report** command again. To cancel the printout request, click on the **Cancel Printout** button.

If a Detailed Report is desired, from either page of the Detailed Results screen, select **File** and click on the **Print Report** command. Then click on either the **Print Immediately** button or **Add to Queue** button. To cancel the printout request, click on the **Cancel Printout** button.

3.2.9 Quantities Measured

When a Summary Results display is requested, values are given for Bone Mineral Density, Bone Mineral Content, and region Length for L2, L3, L4, and the L2-L4 column. The Page 2 Detailed Results display provides the same measurement values given on the Summary Report plus AREA values and CV values ("CV" means Coefficient of Variation, as described in Section 5.6.).

It is generally accepted that in order to have confidence that a change in measurement represents a change in BMD, the measured change must be at least:

$$[(CV \text{ Initial Meas.})^2 + (CV \text{ Final Meas.})^2]^{1/2}$$

⁴

CV values are valid for subjects having well defined bones and non-calcified soft tissue. CV values do not apply to the vertebral segment length figures since these values are determined solely by operator judgement and not XR-26 analysis. Length values are provided for reference only; to assist in repeatable placement of intervertebral cursors on follow-up scans.

⁴LeBlanc, A. D., et. al.: "Precision of Dual Photon Measurements." Journal of Nuclear Medicine. 27:1362, 1986.

3.3 TRENDING

The XR-26 provides several means to visually present a patient's bone mineral density measurement in comparison to several indicators. Please note that **only** the clinician performs the diagnosis. The XR-26 supplies certain data which may assist the clinician.

The Trending technique assists an interpretation of the patient's measurement when compared with previous measurements of the same patient. Trending looks only at the individual patient's measurements. Please refer to Section 5 for a complete description of Trending operations.

The second method of patient comparison is called Reference. All information concerning Reference may be found in the Reference Population supplement to this Operator's Guide.

The descriptions of A-P Spine Scan and Analysis is now complete. We invite you to review other sections of this Operator's Guide.

This page reserved.

4. HIP MEASUREMENT

The hip measurement mode in the XR-26 scans the left or right proximal femur area and provides bone mineral quantification of the femoral neck, Ward's Triangle, and the greater trochanter regions of interest. The analysis software provides a means of excluding scan artifacts and for following a patient serially.

Proper operation of the XR-26 and interpretation of its measurements requires a thorough understanding of the factors affecting precision and accuracy. Such information is provided in Section 5.6, "Precision and Accuracy."

Although a female patient is used for procedure descriptions herein, the same techniques are to be used with a male patient. In addition, references are made to various sections in the Base Operator's Guide. You should become familiar with this referred material before performing the techniques described herein.

The scanner is equipped with a low power helium-neon laser which is used when positioning the scanner arm over the patient. The red laser light spot is projected onto the patient and indicates where the x-ray beam will pass.

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength.
DO NOT STARE INTO BEAM.

4.1 SCANNING THE HIP

In brief, the hip scanning process begins with a quick Scout scan over the femoral neck area. The software creates and displays an image of the scanned anatomy. On this screen image, the operator defines start and end points for the measurement scan. The measurement scan is then taken, results calculated, and the data is saved. During this process, instructions for scanner arm movement and scan site identification is presented as appropriate.

Several differences exist between Scout Scan and Measurement Scan attributes and are intended to provide scan data appropriate to the procedure used. Table IV indicates the differences in fixed scan attributes, as selected at the Verify Scan Parameters screen.

	Scout Scan	Measurement Scan
Point Resolution	1.0 mm	1.0 mm
Line Resolution	2.0 mm	1.0 mm
Scan Speed	80.0 mm/sec	45.0 mm/sec

Table IV - Fixed Hip Scan parameters

4.1.1 Preparing Patient

For a hip scan, the patient need not undress or even loosen clothing. However, there should not be any metal or plastic objects in the scanning area that could interfere with the scan. Check for hip or back braces, any metal on clothing such as rivets, buttons, zippers or other types of fasteners, or metal, coins, or other objects in front or rear pockets. In other words, check for any item that could add a density to the study. If any of these items are found, be certain to keep them clear of the scan area.

The foot block (a trapezoid-shaped fabric-covered foam block) and foot straps should be readily available. This block aids in optimal positioning of the proximal femur for scanning, as seen in Figure 20.

The patient lies down in the center of the patient pad with her right side adjacent to the scanner backrest (this places her head at the 'head' end of the scanner). This orientation is important in hip measurement for the correct identification of left and right hip scans. The scanner arm should be positioned over the patient's midsection using the scanner arm movement keys. Figure 20 shows a typical patient positioning scheme.


Figure 20 - Patient on Scanner with Foot Block

The patient is then arranged so that her legs are placed straight out and her feet are rotated inward slightly. For correct positioning of the feet, place the foot block between the patient's feet and elevate the heels into the block heel rests. Secure the toes to the block using the foot straps, as seen in Figure 21. For your more elderly patients, use the lower heel rests for greater patient comfort. It is good practice to develop a consistent positioning method for all patients.

As the patient may be lying in this position for several minutes, it is important that she is made as comfortable as possible. This care reduces the chances of extraneous movement. Make sure she is warm or cool enough, and that there is a pillow for her head. The use of a sheet or light blanket for warmth or modesty will not interfere with the scan results. The patient is now ready for scanning.


Figure 21 - Proper foot positioning

Select Patient

Please refer to 5.1.1 for detailed instruction on using the patient selection features of the XR-26. If the patient to be scanned is a new patient and needs to be added to the patient directory, please refer to 5.1.1.

4.1.2 Scan Preparation

After selecting the patient, select **Scan**. Notice that separate scan procedures are provided for Left Hip scan and Right Hip scan. Left hip and right hip designation is made from the patient's point of view. These instructions will continue by describing left hip procedures, but are easily applied to right hip procedures by adopting a mirror image of the anatomy presented.

Click on the **Left Hip** command. A screen requesting the patient's height and weight, and the names of the scanning technician and physician appears. This information is entered using Dialog Edit Techniques described in paragraph 1.4.4 and automatically added to the patient's data file. This information is requested each time a patient is scanned.

With the patient in position, the operator should prepare the scanner as prompted by the display. Turn on the laser and observe that a red dot appears directly below the scanner arm. If the red dot is not visible, check that the Laser Manual Shutter paddle is fully elevated (horizontal).

WARNING: The XR-26 is equipped with a Class II Laser which is less than 1 milliwatt in strength. **DO NOT STARE INTO BEAM.**

The typical LEFT and RIGHT proximal femur scan regions are shown in Figure 22. Each region is oriented at a 45 degree angle to the table's long axis and centered about the MARK point. The preset length of a hip scout scan is 4 cm above the MARK point and 6 cm below the scan point. The preset width of the scan is 10 cm.


Figure 22 - Recommended Left and Right Hip Scan Areas

In Figure 22, certain letters are used to identify the scan landmarks. The [M] denotes the recommended Mark point position requested of the

operator. The other landmarks are defined by XR-26 software, where [S] denotes the Start point for the scan, [E] denotes the End point, and [B] denotes the baseline measurement point.

Use a palpable anatomical feature (where the femur is closest to the skin, for example) to locate the proximal femur. Position the red laser dot at the approximate center of the femoral neck area using the arm motion controls. Press the **MARK** key to continue. A beep is emitted by the computer to indicate that the marked point is acknowledged. The **MARK** position doesn't have to be exact, since the ensuing Scout Scan screen image is used to define the Measurement scan. If you click on the **Cancel Scan** button, another message appears, asking if you wish to terminate the operation. Click on the **Resume** button to continue or on the **Terminate** button if you wish to cancel the scan.

Once the center of the hip area has been marked, a Scan Parameters screen similar to Figure 23 appears, listing point resolution, line resolution, scan width, and scan speed parameters. Point and line resolution parameters are preset by XR-26 software and cannot be changed. The Scan Speed is 45.0 mm/sec for Measurement scans and 80.0 mm/sec for Scout Scans.


Figure 23 - Scan Parameters screen

Norland recommends 16.0 cm for Scout Scan Width to provide ample soft tissue characterization. If desired, the scan width may be redefined using the Dialog Edit Techniques, Section 1.4.4. To change the Scout Scan Width, click on the **Scout** button and note that the scan speed value changes to 80.0 mm/sec. Enter the desired scan width. The effect of this change is seen on the Scan Time value displayed on the next screen. This process may be repeated until satisfied with Scan Width and Time for the Scout Scan.

The Scout Scan will be performed using the scan parameters last seen at the Scan Parameters screen, so make sure the proper Scan Speed and Scan Width values are indicated.

If necessary, click on the **Remark Scan** button. Instructions for scanner arm placement appear. Follow the on-screen guidance to re-mark the scan on the patient. Or, if you click on the **Cancel Scan** button, you are asked to confirm your desire to terminate the operation. Click on the **Resume** button to continue or on the **Terminate** button if you wish to cancel the scan.

4.1.3 The Hip Scout Scan

When satisfied with scan parameters, click on the **Continue** button. An X-Ray Parameters screen similar to Figure 24 appears, indicating the intended Anode Voltage (in kilovolts), Anode Current (in milliamps), and Estimated Scan Time (in minutes).


Figure 24 - X-Ray Parameter Information

If you concur with the scan time, caution the patient to lie still, then click on the **Start Scan** button. The XR-26 turns-off the laser and selects the filter combination dictated by patient thickness, then measures a baseline. The computer may display a message stating that the x-ray source is being brought to operating condition. The scanner arm will soon begin to move back and forth across the patient. The system completes the Scout Scan now without intervention, then signals the end of the scan with an audible beep.

If you do not concur with the listed scan time, click on the **Edit Parameters** button. Make the changes you desire to Scan Width or Scan Speed, and click on the **Continue** button. The computer will display the effect of this change on the Scan Time shown in the ensuing X-Ray Parameters screen.

The Scan Progress screen shows how many lines will be scanned and gives an estimate of scan time remaining. While the scanner is acquiring data, an image of the scan is built up on the screen.

It is important for the operator to observe the screen image. It is possible that: the wrong region is being scanned; the scan extends beyond the desired region of interest; the patient moved; the scan image is unrecognizable; or other difficulty may arise. Close observation permits the scan to be terminated, if necessary, thereby minimizing patient radiation exposure.

The image is rendered in a gray scale (or color spectra, if selected). Note the Density Spectrum that appears near the image. The 'L' indicates the lowest detected bone density and the 'H' indicates the highest bone density. Intermediate bone densities are evenly allocated to the intermediate gray (or color) levels. These images aid in placing the cursors and are not intended for diagnostic purposes. Subsequent analysis is performed on the actual data acquired.

Remember that if the **EMERGENCY OFF** key on the Scanner Control Panel is pressed, the scanner must be turned off, then on again. Scanning may then be resumed.

Once scanning has started and the **SCAN PROGRESS** screen is displayed, clicking on the **Cancel Scan** button causes the current scan line to be completed. An information box appears, asking if you wish to terminate the operation. Click on the **Resume** button to continue.


Figure 25 - Scan Complete screen

If you click on the **Terminate** button, a **Scan Terminated** screen appears (very similar to Figure 25). You are offered the same choices as if the scan is completed normally.

4.1.4 The Hip Measurement Scan

When the Scout Scan is successfully completed, the Scan Complete screen (shown in Figure 25) appears. You are asked if you wish to Rescan the patient or to Discard the data. If you click on the Discard button, the data is not saved and the Main Menu - Idle Display eventually reappears.

Click on the **Rescan** button and observe a screen similar to Figure 26. The two "X" cursors connected by a line and a box are superimposed on a reduced-size scan image. Caution the patient to remain lying very quietly. The cursor(s) should be positioned on the image as shown to define the upper and lower limits for the Measurement Scan or another Scout


Figure 26 - Position Rescan Cursors screen

Scan. The choice of Scout or Measurement Scan is made at the Scan Parameters screen, which will shortly be seen. Should you choose another Scout Scan, follow the directions given above.

For measurement repeatability, it is very important to obtain a properly-oriented Measurement scan image, i.e.: the long axis of the femoral neck should be straight up and down on the final measurement scan image. Drag the upper or lower "X" cursor to where you want it on the screen image. Then drag the other "X" cursor to the desired position. The line that connects the two "X" cursors should be used to guide cursor placement so that the ensuing scan is centered on the long axis of the femoral neck. The box outlines the intended scan region.

Once the hip scan cursors have been adjusted, click on the **Edit Scan Parameters** button. A Scan Parameters screen similar to Figure 23 appears. Verify the scan width value shown on the Scan Parameters screen.

Norland recommends 10.0 cm for Measurement Scan Width but it may be re-defined by the operator. To change the Scan Width, click on the **Measurement** button and observe that the Scan Speed is 45.0 mm/sec. Enter the desired scan width using Dialog Editing Techniques (1.4.4).

If necessary, click on the **Remark Scan** button. The software remembers that a Scout Scan has been performed and presents a reduced-size Rescan Cursor positioning screen. Another Scout Scan may be defined there and any work done to this point is discarded. Follow the on-screen guidance to re-mark the scan. Again, make sure the proper Scan Speed and Scan Width values are indicated before leaving this screen. Or, if you click on the **Cancel Scan** button, you are asked to confirm your desire to terminate the operation. Click on the **Resume** button to continue or on the **Terminate** button if you wish to cancel the scan.

When satisfied with scan parameters, click on the **Continue** button. An X-Ray Parameters screen similar to Figure 24 appears, indicating the intended Anode Voltage (in kilovolts), Anode Current (in milliamps), and Estimated Scan Time (in minutes). Review these parameters, caution the patient to lie still, then click on the **Start Scan** button. The scanner arm will begin to move across the patient and the measurement scan will be completed without further intervention. The Scan Progress screen indicates the estimated scan time remaining.


Figure 27 - Measurement Scan Complete screen

When the end of the scan is reached, a beep is emitted and the XR-26 presents the Scan Complete screen; as in Figure 27. You are asked if you wish to Analyze, Rescan, or Discard the data. Verify the proper orientation of the scan image.

If you terminated the scan, a Scan Terminated screen very similar in appearance to Figure 27 is displayed. You are offered the same operational choices as if the scan was completed normally. If you click on the **Rescan** button, you may perform a rescan of the patient, as described above. If you click on the **Discard** button, the data is not saved and the Main Menu / Idle Display reappears.

Move the scanner arm out of the way and detach the foot block from the patient. The patient may now be released.

4.2 HIP ANALYSIS

The information presented here continues with the standard scan-then-analyze sequence of operations. XR-26 software provides a way to re-analyze scan data and is described later in this section.

In a Hip Analysis, the characteristics of the femoral neck, Ward's Triangle, and the greater trochanter are calculated based on the information stored in the scan data file. The hip analysis process begins with cursor positioning to define regions. If scan artifacts need to be removed from analysis consideration, an Include / Exclude process is provided. Numeric results are then calculated and presented.

To proceed to Analysis, click on the **Analyze** button.

4.2.1 Identify Femoral Neck

The analysis process begins with placement of the femoral neck cursor so that it cuts the left window neck image at its narrowest point. Cursor placement during this process is largely automatic, but at each step, the operator is asked to approve the program's choice and is allowed to override it if desired.

The Left Hip Cursor Placement display (Figure 28) shows the various bone densities found during the scan in a grayscale or color density spectrum. 12-level gray scale, 16-level gray scale, or 12-level color scale image depiction may be made by selecting **Image** and clicking on the desired rendition command.

For a routine hip measurement, the scanned region should include the entire femoral neck, the entire greater trochanter, and part of the

femoral shaft and the ilium. The center line of the scan path should coincide with the center line of the femoral neck.

It may prove useful to imagine a straight line that bisects an image of the femoral neck and is parallel with the neck. Then imagine another straight line that is perpendicular to the first line and passes through the center of the femoral neck. The closer the initial MARK point is made to the intersection of these two lines, the better the final measurement will be. Consistent use of this visualization when marking the screen image should result in successful and repeatable analyses.


Figure 28 - Femoral Neck Cursor Positioning screen

The XR-26 software analyzes the data and automatically places the horizontal neck cursor to contain the narrowest segment of the femoral neck. The operator is asked to verify this positioning and, if necessary, change it.

If you do not concur with automatic femoral neck cursor placement, move the pointer near the center of the neck cursor and drag the entire cursor to the narrowest portion of the femoral neck.

A single vertical line is positioned at the center of the neck segment defined by the bone edges and horizontal cursors. If this line appears entirely outside the bone edges of the femoral neck, drag it into the neck segment between the bone edges.

The XR-26 aids the operator by moving the cursors together as a group; by displaying the neck length on the screen (and on printouts) for reference; and by offering a Comparison Image. Improper cursor positioning may render the analysis results inappropriate for serial studies.


Figure 29 shows an example side-by-side Comparison Image. To obtain a similar Comparison image, select *Image* and click on the *Show Comparison* command. The analyzed image from the first left hip scan for the patient is recalled and presented in the center of the display, using the same linear scale as the existing image.

Note: When performing follow-up scans on patients, it is very important that the cursor placement and size be the same as the initial scan. In particular, the neck cursor height should be the same. This dimension is displayed on the screen for reference. The operator is encouraged to use the Comparison Image or to refer to the printout of the initial scan when verifying the cursor positions.

Baseline data (non-bone points) are shown or 'hidden' by clicking on the appropriate *Image* command. Removal of baseline data reveals all scanned points corresponding to the chosen display spectrum without regard to their 'bone' or 'non-bone' status. This feature may be helpful in locating anatomical landmarks such as ribs and iliac crests, which have been automatically excluded from the analysis. To select an alternate image rendition, select *Image* and click on the rendition command of your choice.

If artifacts have been included as a part of the analysis, they may be removed using the INCLUDE / EXCLUDE feature, as described in Section 4.2.4.

4.2.2 Numeric Results

When satisfied with both cursor positions, select **Results** and click on the **Summary Results** command. The program locates the Ward's Triangle marker, performs summary calculations, and stores the results. A tabular report of the bone measurements is then presented, as shown in Figure 30. The image is not intended for diagnostic purposes.


Figure 30 - Summary Results screen

The Summary Results screen includes the scan image with femoral neck and Ward's Triangle cursors in position. Listed are measurement results for bone mineral density (BMD), bone mineral content (BMC), and neck segment length (LGTH). An estimate of the coefficient of variation (CV) for each measurement is included in parentheses beneath it.

XR-26 Software provides two additional Results screens. The Detailed Results page 1 screen, seen in Figure 31, presents the scan image alongside the patient personal data. The Detailed Results page 2 screen, seen in Figure 32, lists the values measured for each proximal femur region and associated CV values.


Figure 31 - Detailed Results Screen, Page 1



Figure 32 - Detailed Results Screen, Page 2

Comments about this scan, such as special conditions or observations, may be entered using Dialog Edit Techniques (1.4.4) after selecting **File** and clicking on the **Edit Comments** command. Any comment so entered will be saved with the data and is included on any subsequent numeric results printout. Click on the **Accept Changes** button to preserve the scan comments.

4.2.3 Reviewing Previous Results

When an analysis for a particular patient is desired, it is necessary to specify the patient and the scan of interest. To do so, select **File** and click on the **Select a Patient** command. Double-click on the particular Patient name, then double-click on the desired Left or Right Hip scan data file. The computer then retrieves the data for the selected scan and presents a Results screen. More detail about the patient or scan selection features of XR-26 software may be found at 5.1.1.

If the data was stored on diskette, a "Loading Data" message requesting the insertion of that specific Volume disk is displayed. Place the diskette in the disk drive and click on the **Continue** button.

For Summary and page 1 Detailed Results, neck cursor positions are recalled with the data and are displayed in position. If cursor positioning needs to be redefined, select *Regions* and click on the *Modify Regions* command. The Review Cursor Placement display is presented and the desired changes may be made, as described in paragraph 4.2.1. If additional artifacts have been identified, select *Regions* and click on the *Include/Exclude* command. Mark and process the artifacts as described below.

4.2.4 Artifact Include / Exclude

The hip analysis program has been designed to automatically exclude any 'bone points' which are not connected to the pelvis or femur. Similarly, it will automatically eliminate any voids or islands of soft tissue surrounded by bone. It may be sometimes necessary to manually remove non-femur artifacts from the measurement region by using the Exclude function. To do so, from the displayed Results screen, select *Regions* and click on the *Modify Regions* command. Reselect *Regions* and click on the *Include/Exclude* command.

On the Include/Exclude screen, a small box appears over the image (as in Figure 33). It will be used to circumscribe such artifacts for removal from analysis consideration.

Position the pointer anywhere inside the box and drag the box to the position desired. Then, drag an edge of the box where desired. Continue with this positioning sequence for the remaining sides of the cursor. When satisfied with the


Figure 33 - Include/Exclude screen

box size and placement, click on the **Exclude** button to exclude the region from results consideration.

If you wish to include a portion or all of an excluded region, click on and drag each edge of the cursor box to circumscribe the area of interest. When satisfied with cursor size and placement, click on the **Include** button. Data from the newly-defined region is included in any subsequent results report.

This process of excluding each fragment of the artifact may be repeated as often as needed. When satisfied with region exclusion, click on the **Continue** button. Review the neck and center cursor positions, then select **Results** and click on the desired results command. Data from any excluded region(s) is removed from the tabulated results. Data from any included region(s) is added to the tabulated results.

4.2.5 Rescale Image and Reanalyze Scan

It is possible that a scanned region includes a very high density object. The Report image density spectrum is scaled according to the range of densities found in the scan. An high-density object included in the analysis causes the resultant density spectra to be spread across a wider-than-normal range of density values, so the on-screen density rendition may be skewed. The portrayal of lower-density scan elements may be obscured because of this skewed density spectrum.

XR-26 software permits these high-density objects to be removed from analysis consideration using the Exclude feature. Once removed, the scan image may be rescaled so that a more normal range of density values is portrayed. To do so, use the instructions given above to exclude the object. Then select **Image** and click on the **Rescale Image** command. The software will disregard the excluded areas from the analysis and present the image in a more realistic density spectrum.

The Include feature may be used to recover each excluded portion of a scan in a piece-by-piece fashion. The Reanalyze operation may be used to recover all excluded regions in one operation, without moving

any Hip or Special Region cursor(s). To reanalyze a scan and recover all excluded regions, select *Regions* and click on the *Reanalyze* command. Confirmation of this operation is NOT made, so be sure you have no need for any of the existing excluded regions. The analysis process starts over from the beginning.

A High Density Detection Option is available for XR-26 software. If installed and enabled, the software will discard any data point containing a very high density value, thereby removing these artifacts from the computed measurement values. This Detection Option may be enabled or disabled at the Analysis Results screen by selecting *Regions* and clicking on the desired *High Density Detect* command. Once set, this mode remains in effect until changed.

4.2.6 Printed Report


Figure 34 - Left Hip Analysis Printed Report

Figure 34 is an example of a typical Left Hip Summary printed report. The printed image is not intended for diagnostic purposes.

Choices made in the *Image* Menu are used when generating hard-copy printout of the analysis. These choices include images rendered with or without baseline data and in color or gray scale. What is seen on the display is what will be printed. The number of copies that will be printed is pre-selected at Printer Setup, described in Section 5.4.4.

Two types of printed reports may be obtained. The Summary Report, shown in Figure 34, depicts essentially the same information seen at the Summary Results screen. It takes less time to generate and print than the Detailed Report, which combines the information seen in both Detailed Results screens on one piece of paper.

If a Summary report is desired, from the Summary Results screen, select **File** and click on the **Print Report** command. Then click on either the **Print Immediately** button or **Add to Queue** button. The preselected number and rendition of the analysis results is obtained. If additional printed reports are desired, click on the **Print Report** command again. To cancel the printout request, click on the **Cancel Printout** button.

If a Detailed Report is desired, from either page of the Detailed Results screen, select **File** and click on the **Print Report** command. Then click on either the **Print Immediately** button or **Add to Queue** button. To cancel the printout request, click on the **Cancel Printout** button.

4.2.7 Quantities Measured

When a Summary Results display is requested, values are given for Bone Mineral Density and Bone Mineral Content for the Neck, Trochanter, and Ward's Triangle regions, as well as the Neck Length. The Page 2 Detailed Results display provides the same measurement values given on the Summary Report plus neck, trochanter, and Ward's Triangle AREA values.

It is generally accepted that in order to have confidence that a change in measurement represents a change in Bone Mineral Density, the measured change must be at least:

$$[(CV \text{ Initial Meas.})^2 + (CV \text{ Final Meas.})^2]^{1/2}$$

5

"CV" (Coefficient of Variation, as described in Section 5.6.) values are valid for subjects having well defined bones and non-calcified soft tissue. CV values do not apply to special region figures since these values are determined solely by operator judgement and not XR-26 analysis. These values are provided for reference only; to assist in repeatable placement of cursors on follow-up scans.

⁵LeBlanc, A. D., et. al.: "Precision of Dual Photon Measurements." Journal of Nuclear Medicine. 27:1362, 1986.

4.3 TRENDING

The XR-26 provides several means to visually present a patient's bone mineral density measurement in comparison to several indicators. Please note that **only** the clinician performs the diagnosis. The XR-26 supplies certain data which may assist the clinician.

The Trending technique assists an interpretation of the patient's measurement when compared with previous measurements of the same patient. Trending looks only at the individual patient's measurements. Please refer to Section 5 for a complete description of Trending operations.

The second method of patient comparison is called Reference. All information concerning Reference may be found in the Reference supplement to this Operator's Guide.

The descriptions of Left and Right Hip Scan and Analysis procedures are now complete. We invite you to review other sections of this Operator's Guide.

5. GENERAL TOPICS

5.1 SPECIFYING SCAN/ANALYSIS SUBJECT

When the XR-26 program begins its operation, it always pre-selects the first scan of the first patient found in the patient directory. The operator may use the following techniques to select another patient or scan before performing new scans or analyses.

5.1.1 SELECT A PATIENT

When **File** is selected from the Main Menu, a drop-down menu appears listing various patient record operations. Click on the **Select a Patient** command, and the SELECT A PATIENT screen appears. Figure 35 is a sample listing of patients in the patient directory.


Figure 35 - Sample Patient List screen

To select a patient from this listing, simply move the pointer until it is on top of the patient's name. To display the first screen-full (or last) of patient names, click on the **Home** (or **End**) button. Clicking on the **PgUp** (or **PgDn**) button displays the previous (or next) screen-full of names.

Two search functions are provided and may be used to help find a particular patient. To search for a patient by name, simply click on the **Search by Name** button, enter the first few characters of the patient's name, then click on the **Begin Search** button. Use either upper or lower case letters. The computer will search the entire patient list for the first occurrence of a name that matches the characters you entered.

Depending on the number of characters used to define the search string, several similar names may be encountered before the desired patient record is accessed. If no match is found, the patient who's name most closely matches the Search string appears in the center of the patient list screen.

The search of the patient list for a particular I.D. is very similar to the method described above. Click on the **Search by ID** button, enter the first few characters of the desired patient's ID number, then click on the **Begin Search** button. The computer will search the list for the first ID number that matches the characters entered. The patient record thus found will appear in the center of the patient list screen.

Leaving this display causes the highlighted name to be 'selected'. This means the name and ID is displayed at the top of the screen for all subsequent screens until another name is chosen.

Adding New Patient

When the patient to be scanned is a new patient, select **File** and click on the **Add a Patient** command. The Add a Patient display appears (Figure 36), requesting that the name and I. D. number of the patient be entered using the Dialog Edit Techniques (1.4.4). The last name may be recorded first, then a space, then the first name or initial.

When the name entry is satisfactory, press the **Tab** key (not the **Enter** key) and the cursor moves to the I.D. line. Use the keyboard to enter the identification number desired. Norland suggests that some sort of unique number be used, like the patient's hospital registration number, birth date, or a specially assigned clinic patient number. In any case, it

is necessary to enter something in each box. Note that the Tab and Shifted-Tab keys move the cursor back and forth between the boxes.

If either the name or the I.D. so entered is the same as an existing patient on the patient list, a message appears saying that this patient already exists. Re-enter a different name or I.D. number, or click on the **Do Not Accept** button to return to the Select A Patient screen.

Figure 36 - Add A Patient screen

Once the patient information has been successfully entered, the Select A Patient list appears with the new name highlighted and selected (in this case, selected means the name/ID appears at the top of the display and is used by the system until another name is chosen).

Patient Personal File

Each patient in the patient directory includes a personal data file. This record (Figure 37) contains entry fields for such items as address, telephone number, birth date, bone history, treatment, billing information, and other patient record items. Certain fields, like Menopause Age, must contain integers. Attempts to enter characters in these fields will be ignored. This information appears on

Figure 37 - Patient Personal Data edit screen

Analysis Results printouts and the page 1 Numeric Results screen, and may be updated at will.

Use the mouse and Dialog Edit Techniques (1.4.4) to enter the desired information. The **Tab** key moves the entry cursor to the next field, shift-tab moves the cursor back to the previous field. When all information is as desired, click on the **Accept Changes** button. If you click on the **Do Not Accept** button, the information previously entered is discarded. Pressing the Enter key is the same as clicking on the **Accept Changes** button.

The Birthdate field must contain information if any kind of Trending or Reference plot is to be made. Depending on the preset form of date, the birth date is entered for MM=month, DD=day, YY=year. Should the date be missing when a Trending chart is requested, the scan dates are used instead of Birthdate.

Delete a Patient

If it becomes necessary to remove a patient from the patient directory, the patient records should be copied onto archive diskettes using the copy process described in Section 5.3.3. Unless they are copied, patient records are not recoverable after the deletion process.

To delete a patient record, click on the name of the patient in question, then select **File** and click on the **Delete a Patient** command. The Delete Patient confirmation box (Figure 38) appears. Clicking on the **Delete** button causes all scan data and personal records for the selected patient to be erased, including that data stored on diskette. When appropriate, screen messages request that certain Volume diskettes be inserted for the deletion process. Clicking on the **Cancel** button causes the process to stop without affecting the patient data and the Idle screen appears.


Figure 38 - Delete Patient confirmation box

5.1.2 SELECTING SCAN DATA FILE

The list presented at the Select A Scan display is a record of all scans taken for the patient in question. If this patient had been scanned three times, each scan would be listed here. Any one of the scans may be selected for an operation. Each scan entry indicates the type of scan taken, when it was taken, and the Volume Number of the disk where the scan data was stored. Data is stored on the hard disk if the Volume Number is '0'.

If more than one scan was performed for the patient on the same date, the sequence number listed next to the scan date is incremented. Thus, a second scan taken on any particular date would have a '2', a third scan a '3', and so forth.

To select a particular scan for re-analysis, double-click on the desired scan. The next screen to appear depends on which Results display was last used, and further re-analysis operations may be performed.


To delete a scan record, click on the scan in question, then select **File** and click on the **Delete a Scan** command. The Delete Scan confirmation box asks that this deletion process continue.

5.2 SCANNER

The operator may initiate the daily Calibration procedure, direct the scanner to find the table dimensions (Active Scanning Area) or the table origin, or review the previous XR-26 instrument calibration using the commands that appear in the **Scanner** drop-down menu. The **Begin Calibration** command is thoroughly discussed in Section 2.3.1. The other operations are described below.

5.2.1 REVIEW CALIBRATION

Clicking on the **Review Calibration** command presents the results of the latest instrument calibration. A sample Review Calibration Report screen is seen in Figure 39. This data is derived from measurements made on the Calibration Standard supplied with your XR-26 system.

The screenshot shows a window titled "XR-26" with a menu bar containing "File", "Scan", "Analysis", "Scanner", "Disk", "Utilities", and "Setup". Below the menu bar, it says "Patient: Test Scans" and "Scan: General 4 on 06/28/98". The main content area is titled "CALIBRATION RESULTS" and contains the following data:

		Polynomial	Residuals	
			Metal	Plastic
Last Calibration:		Filter 0	0.0036	0.0074
06/29/98 14:56		Filter 1	0.0036	0.0071
Software Revision: 2.1.0		Filter 2	0.0034	0.0065
Scanner Revision: 1.1.4		Filter 3	0.0023	0.0035
Standard ID: 001		Filter 4	0.0026	0.0046
K1: 0.43946		Filter 5	0.0028	0.0050
K2: -0.00262		Filter 6	0.0029	0.0049
K3: -0.02642				
K4: 3.03752				
K5: 1.53572		Rough Body	0.0510	0.0776
		Spline 0	0.0029	0.0049
		Spline 1	0.0121	0.0200
		Spline 2	0.0213	0.0360
		Spline 3	0.0125	0.0203
		Spline 4	0.0151	0.0232
		Spline 5	0.0401	0.0607
Discriminators: 120 70				
Background: 1 2				

Figure 39 - Review Calibration Report screen

If a printed Calibration report is desired, select **File** and click on the **Print Report** command. Then click on either the **Print Immediately** button or the **Add to Queue** button. The requested report is obtained. If additional printed reports are desired, click on the **Print Report** command again. To cancel the printout request, click on the **Cancel Printout** button.

5.2.2 FIND TABLE ORIGIN

During routine operation of the XR-26, it is possible that the scanner arm may be bumped or nudged hard enough to cause the mechanical drive belts to 'slip a notch' or otherwise be dislodged. Should sufficient dislocation occur, the software may 'lose track' of where the arm is physically located and may cause an error message to appear. In this event, or at any time the operator feels that the arm has been jarred or dislodged, the Find Origin command should be used.

To initiate a Find Origin command, select **Scanner** and click on the **Find Table Origin** command. The scanner arm will move to the foot end of the table, gently easing into the rear corner, to determine the physical home position. When complete, this origin value is stored as the HOME value and remains in effect until another request to find origin is made.

5.2.3 FIND TABLE DIMENSIONS

Only a small portion of the perimeter of the scanner pedestal is used for mechanical support of the scanning surface, thus providing a large active scan area. When the instrument is first installed, the Norland Installer will instruct the scanner to determine the exact size of the active scanning area using the Find Table Dimensions command. The values acquired with this operation are stored and used to verify that a requested patient scan procedure will fit entirely within the active region.

To establish the scan limits, select **Scanner** and click on the **Find Table Dimensions** command. The scanner arm will traverse the length and width of the table to determine the physical limits of arm travel. When complete, the acquired values remain in effect until another request for table limits is made.

5.3 DISK UTILITIES

The XR-26 and XR-26 Mark II systems now permit the use of multiple hard drives, multiple floppy drives, optical disks, and other DOS-compatible data storage devices, for storage of XR-26-acquired patient data. These capabilities come with the purchase of the Configurable Disk Drives option. The operation and use of this option is fully explained in the Configurable Disk Drives Supplement.

To simplify data storage explanations throughout this Operator's Guide, all removable media types will be referred to as 'diskettes'. And all hard drive operations will be described as 'C:' drive operations, even though you may use other drive designators for your scanning operations.

Disk Memory Philosophy

Rule Number 1. The XR-26 system is designed so that the program itself and all patient scan lists are preserved on the hard disk. The more lengthy and less frequently accessed scan data files are kept on Volume diskettes.

Rule Number 2. Since the patient scan records are irreplaceable and not recoverable, Norland strongly recommends that routine hard disk backup procedures be instituted to store the most recent patient records.

Although there are NO exceptions to Rule Number 2, there are exceptions to Rule Number 1. One is that the operator may choose to keep some scan data on the hard disk for ready reference, with the understanding that these files rapidly use up hard disk space. Another

exception may be if an on-going research task has dedicated the instrument to that task.

In these cases, any one of the popular Disk Operating System hard disk backup utility programs permits a rapid hard disk backup procedure. The only limitation of this procedure is the potential loss of scan data acquired since the last time the backup procedure was performed. Obviously, the regular use of a backup procedure eliminates the need to perform the backup procedures described elsewhere in this operation guide.

Storing Scan Data

While a scan is in-process, the acquired count data is converted to its bone and soft tissue components. It is stored at the end of the calculation process in disk memory. This disk memory is either the internal hard disk drive or external diskettes, as selected by the operator. This choice is made from the Setup menu, as described in Section 5.4.3. External diskettes which are part of the XR-26's data storage system are called "Volume" diskettes because each of them has been assigned a volume number.

If the proper volume diskette is in place, the scan data is written to it without operator intervention. If REMOVABLE Media has been selected as the scan data storage area and a Volume diskette is not present in the disk drive at the time data is ready to be stored, the software will present a message asking if the data is to be stored on the hard disk (FIXED Media). Click on the OK button to initiate storage to the hard disk, and the process is complete. Or, if the data is to be stored on a volume diskette, insert the proper volume diskette before proceeding.

If a new formatted diskette is inserted, the computer will record on it the next available volume number before storing the scan data. This process is called 'disk initialization'. If a proper volume diskette is in place, but does not have enough space left on it to hold the scan data file, a message appears asking that a different diskette be inserted.

5.3.1 DISKETTE STATUS

It is possible that the written Volume diskette labels may become damaged or separated from the diskette. The XR-26 program permits the inspection of a diskette for a Volume label file by pressing a few keys. To do so, insert the diskette in question into the disk drive. Select *Disk Utilities*, click on the *Diskette Status* command, and observe the status message that appears on the display.


Figure 40 - Diskette Status information box

An example of the Diskette Status message for a Volume Diskette is shown in Figure 40. Notice that 265 kilobytes is the remaining disk storage space on this Volume 01 Diskette. Similar messages are shown for the other types of media present in the requested drive.

5.3.2 MOVING DISK FILES

Operations such as Delete Scan or Patient, Move Scans, and Copy Scan or Patient, have been provided for handling the manipulation of patient disk files. Except for Delete Patient, described earlier, these commands are described in the order they are listed on the Disk Utilities drop-down menu.

Move A Scan

In XR-26 Software, "Move" has the specific meaning of transferring a specified scan file to another disk storage location. A patient's scan data file may be moved from the hard disk to a Volume diskette or from

a Volume diskette to the hard disk. This function does not make a copy of a file, but relocates the file from one storage location to another.

The process for moving a file is straightforward. Use the Select a Patient routine to identify the patient and scan record you wish to move. Then select **Disk Utilities** and click on the **Move a Scan** command.

An error message is displayed if a requested destination Volume diskette is not present in the removable media drive. After the disk is inserted, click on the **Retry** button to initiate the Move process. 'Move' copies the specified scan file to the destination storage media, then removes the scan data from the original storage media.

It is not necessary to use a new Volume diskette for the 'Move' operation. Any Volume diskette may be used, provided there is enough room remaining for the patient files being moved. The program will check the remaining space on Volume diskettes before trying to move the information.

If a newly formatted diskette is inserted, the XR-26 asks if a volume number is to be written to it. If an incorrect diskette has been inserted, this is the time to replace it with the appropriate diskette. Otherwise, click on the **Continue** button. The diskette is then labeled with the next available volume number and the requested file transfer takes place. When the move operation is complete, the Idle Menu reappears.

Move All Scans

This Disk Utilities menu option permits convenient on-demand transfer of patient scan data from the hard disk to Volume diskettes. Some clinical sites may choose to store all patient scan data on the hard disk during the day's operations in order to save time. Then, at the end of the day, the operator may transfer all the scans to volume diskettes using the Move All Scans procedure. This process will remove scan data from the hard disk after the files have been recorded on the Volume diskette.

To move all scan files, select **Disk Utilities** and click on the **Move All Scans** command. Once initiated, this unattended task simply requests that formatted media be inserted in the media drive. Click on the **Continue** button when ready. The diskette is initialized into a Volume diskette, if necessary, then all patient scan data stored on the hard disk is transferred to the Volume diskette. Should there be more scans to move than can fit onto the Volume disk, the operation will pause and request that another Volume disk be inserted.

The Volume diskettes used by the XR-26 system are 'machine specific'. This means that the data is safely accessible only to the instrument that created the Volume diskette itself. Copy diskettes were created to safely transfer patient data from one instrument to another.

5.3.3 COPYING DATA FILES

The XR-26 **Copy** command has the specific meaning of copying a scan data file, plus all of the associated patient personal history information, onto a floppy diskette. **Copy** is normally used to archive inactive patient files or for sending patient scan data to another XR-26 site.

The copy process makes a duplicate record of the specified scan file from the Volume diskette (if that's where the data currently resides) onto the Copy Diskette. The original Volume diskette or hard disk data is not affected by the copy process.

The Copy diskettes used in this procedure do not have volume numbers associated with them. In fact, the XR-26 will not copy patient scan data onto a Volume diskette. It is possible to have many patients on a single copy diskette (due to the varying size of scan records.)

To copy scan files, use the Select a Patient routine to specify the patient whose file(s) you wish copied. Then select **Disk Utilities** and click on the **Copy a Patient** command. The process is automatic and will display appropriate messages for exchanging diskettes.

If only one or two files for the selected patient need to be copied, use the **Select a Scan** routine to specify the file of interest. Then select **Disk Utilities** and click on the **Copy a Scan** command. Screen messages will appear requesting specific operator actions. Follow the instructions to complete the task.

When the Copy process is complete, the Disk Utilities menu reappears. Repeat the above procedure to copy additional scan records for the same patient or for other patients.

5.3.4 GET COPIED FILES

This command is used to retrieve files originally copied onto Copy or Volume diskettes. Three ways to retrieve scan data are provided.

Get Copied Scan

If only one or two scans need to be retrieved, select **Disk Utilities** and click on the **Get Copied Scan** command. The program requests that the Copy diskette containing the desired scan files be inserted in the disk drive. The program reads the disk and lists the patients and scans that are contained on the disk. Click on the scan you wish to retrieve, and the patient name and I.D. is cross-checked against all the patients already contained in the XR-26 system's Patient Directory.

If no match is found, a new patient sub-directory is created on the hard drive, and the requested patient scan file is copied. The operator may Move the scan files from the hard disk onto volume diskettes at a later time.

If the Patient Directory contains a patient whose name and I.D. matches the selected record, the subject files are copied on the proper Volume diskette. No attempt is made to identify individual scans, so it is possible that duplicate patient scans may exist after the copy process. Repeat the above procedure to copy other patient records.

Get Copied Patient

To retrieve all scans for a single patient, select **Disk Utilities** and click on the **Get Copied Patient** command. The screen requests that the Copy diskette containing the desired files be inserted in the disk drive. From the ensuing patient listing, click on the desired patient. The patient's name and I.D. is cross-checked against all the patients already contained in the system's Patient Directory.

In a manner similar to Get Copied Scan, the requested patient scan file is copied to the hard disk. The operator may Move the scan files from the hard disk onto volume diskettes at a later time.

Get All Copied Scans

To retrieve all scans residing on a Copy diskette, select **Disk Utilities** and click on the **Get All Copied Scans** command. The screen requests that the diskette containing the desired files be inserted in the disk drive. The Copy program then automatically cross-checks the patient list against all the patients already contained in the system's Patient Directory. It creates new directories when needed and adds all files to the hard drive. The operator may Move the scan files from the hard disk onto volume diskettes at a later time. When the Get All Copied Scans operation is complete, the Idle Menu reappears.

The only difference between Volume and Copy diskettes is that Volume diskettes have a volume number file written to them. Copy diskettes do not.

5.3.5 DOS OPERATIONS WITHIN XR-26

Descriptions of XR-26 operations for preparing removable media follows. Complete instructions for the use of similar DOS commands are provided in the Disk Operating System manual.

Format a Diskette

The operation called 'formatting' prepares unused or previously-used storage media for use as disk memory. This format process destroys all information that may be contained on a diskette. We recommend that the Diskette Status command (in the Disk Utilities menu) be used to check previously-used diskettes to see if valuable information could be destroyed.

To format a blank diskette for use as a copy or volume disk, select **Disk Utilities** and click on the **Format a Diskette** command. The XR-26 program asks you to verify the Format process by clicking on the **Continue** button. The screen clears and a message asks that you insert a diskette in the preselected drive. Follow the displayed instructions to complete the process. You may format as many diskettes as you wish. Your "n" response to the "Format another diskette" message will return you to the XR-26 program. If you wish to cancel the format operation without formatting any diskettes, hold down the **Ctrl** key and tap the **C** key. Answer the resulting DOS message as desired.

Backup Files

The XR-26 patient and scan files are extremely difficult to reconstruct in the event of a hard disk media failure. Therefore it is in the operator's best interest to develop and adhere to an effective disk memory backup procedure. Whenever XR-26 software is operating, a hard disk backup of these important system files is maintained. An easy to use XR-26 system file backup and restore procedure is described here.

To backup the important Setup, Calibration, and other system files, select **Disk Utilities** and click on the **Backup Files** command. A "Backup files Continue Cancel" message box appears. Click on the **Continue** button to proceed. The computer screen will display further instructions.

If you wish to cancel the Backup operation, hold down the **Ctrl** key and tap the **C** key. Answer the resulting DOS message as desired. You will be returned to the XR-26 program.

Otherwise, follow the on-screen instructions, and all of the important XR-26 files are backed-up on the diskette. If any errors are detected in the backup process, messages will be displayed and should be reported to Norland Customer Service.

Copy a Diskette

During the sequence of XR-26 operations, it may become necessary to share Copy or Volume diskettes with others. This XR-26 command permits a complete disk copy operation identical to the DOS DISKCOPY program. Please refer to the DOS manual for more information on DISKCOPY.

To make duplicate copies of diskettes, select **Disk Utilities** and click on the **Copy a Diskette** command. A "Copy a Diskette Continue Cancel" message box appears. Click on the **Continue** button to proceed. The computer screen will display further instructions.

If you wish to cancel the disk copy operation, hold down the **Ctrl** key and tap the **C** key. Answer the resulting DOS message as desired. You will be returned to the XR-26 program.

Otherwise, follow the on-screen instructions, swapping diskettes as requested. All the files on the Source diskette are copied to the Target diskette. If any errors are detected in the disk copy process, messages will be displayed and should be reported to Norland Customer Service.

Verify Diskette Copy

The Verify Diskette Copy process is identical to the DOS DISKCOMP program. Please refer to the DOS manual for further information about DISKCOMP. This XR-26 operation will direct the swapping of two diskettes in the preselected drive, comparing the two diskettes section by section. Should any differences be detected, appropriate messages will be displayed.

To verify that two diskettes have identical contents, select **Disk Utilities** and click on the **Verify Diskette Copy** command. A "Verify Diskette Copy Continue Cancel" message box appears. Click on the **Continue** button to proceed. The computer screen will display further instructions.

If you wish to cancel the disk compare operation, hold down the **Ctrl** key and tap the **C** key. Answer the resulting DOS message as desired. You will be returned to the XR-26 program.

Otherwise, follow the on-screen instructions, swapping diskettes as requested. All the files on the First diskette are compared to the files on the Second diskette. If any differences are detected between the two diskettes, messages will be displayed and should be acted upon as desired.

5.3.6 RESTORE HARD DISK

At most locations, all program and data files are preserved on diskettes, including the patient records and other XR-26 files. Each time XR-26 is started, vital system files are copied onto separate parts of the hard disk. These same system files are copied to the BACKUP floppy disk each time the Backup Files command is executed. These files will be used to restore the XR-26 system should that information be lost or destroyed.

Before restoring the information to a repaired hard disk drive, contact the Norland Customer Service technician for detailed instructions. It is a simple process to restore the computer to full working order. Be advised that the backup procedures described here safeguards only XR-26 programs and data. Other programs and data which may reside on the computer's hard disk are not saved. Norland recommends that a commercially-available back up program be purchased and used if other programs or data are kept on the hard disk.

This completes the explanations about the standard disk and diskette use in the XR-26 system. If you have purchased the Configurable Disk Drives option, then descriptions of additional disk operations may be found in the Configurable Disk Drives Supplement to this Operator's Guide. If any further assistance is needed, please review your DOS manual or contact your local Norland representative.

5.4 SETUP OPERATIONS

All operator-selected preset conditions are found under the Main Menu bar title **Setup**. This menu provides access to the initial setup operations for essentially all input and output processes of the XR-26 instrument.

Select **Setup** from the Main Menu bar titles and observe the appearance of a drop-down menu similar to Figure 41. These commands are described below in the order that they appear on the menu. To select the desired setup operation, click on its command. To return to the Idle screen, simply click on any open area of the screen.


Figure 41 -Setup Command Drop-down Menu

5.4.1 CALIBRATION

Clicking on the **Calibration** command presents a display similar to Figure 42. The Calibration Standard Identification values are normally entered by the Norland Installer and should only require re-entry if the standard is changed or if the Setup file becomes corrupted. The Serial (Reference) Number and 'K' values may be found on a label affixed to the Calibration Standard.

When this screen appears, enter the calibration standard's Reference number as the Serial Number. Press the **Tab** or **->** cursor key or the pointer to move the entry cursor to the 'K1' field. Type in the K1 value using the numeral keys, then press **Tab** to move to the next field.

Continue this process to enter the remaining 'K' values.

Once entered, these values remain in effect until changes are made by the operator. This completes the Calibration Setup. Click on the **Accept Changes** button to return to the Setup menu or click on the **Do Not Accept** button to discard the numbers you just entered.

CALIBRATION STANDARD IDENTIFICATION

Serial Number: _____

Calibration Factors

K1: 0.43946_ K2: -0.00262_
K3: -0.02642_ K4: 3.03752_
K5: 1.53572_

Figure 42 - Calibration Standard data entry screen

5.4.2 DATE FORMAT

The operator may select the format (the way that it appears) that the date is displayed and printed during XR-26 operations. The date may be in either the 'MM/DD/YY' form (where MM is month, DD is day, and YY is the last two digits of the year), or the 'DD/MM/YY' form.

The Date Format screen (Figure 43) is where the form of date is selected. To access this screen, select **Setup** and click on the **Date Format** command. When the information box appears, click on the form of date you desire, then click on the **Accept Changes** button. The next time the date is seen, it will appear in the form selected. If you click on the **Do Not Accept** button, the previously-set Form of Date is retained, and the Setup menu reappears.

FORM OF DATE

MM - Month
DD - Day of Month
YY - Year

Figure 43 - Form of Date Selection screen

Setting Date and Time

One important aspect of collecting patient data is to keep track of when that data was acquired. XR-26 Software automatically takes care of this detail by routinely applying a date and time stamp to any disk file it stores. The date and time information is obtained from a non-volatile part of the host computer's memory. Even when the computer is turned off, the date and time data is preserved. Therefore, it is important to make sure this information is correct.

Date and time is set for XR-26 software directly from the computer's 'C>' prompt. Simply type the word 'DATE' and press the **Enter** key. The screen shows the date as it exists in the computer's memory. To reset the date, enter two digits for MM (month), the (/) slash key, two digits for DD (day), the (/) slash key, and two digits for YY (last two digits of the year). Press the **Enter** key to complete the date entry.

To reset the time, simply type the word 'TIME' and press the **Enter** key. The screen shows the time as it exists in the host computer's memory. The desired time is entered as HH:MM:SS, where HH is hour (twenty-four hour clock), MM is minute, and SS is second. Include the (:) colon separator as shown. When satisfied with the time entry, press the **Enter** key.

5.4.3 DEFAULT STORAGE

The XR-26 program provides choices of where acquired scan and results data is to be stored. As seen in Figure 44, the Standard software permits a Fixed Media choice of the 'C:' drive and a Removable Media choice of the 'A:' drive. Note that both choices are highlighted. As described earlier, a Configurable Disk Drive option is available for XR-26 systems that permits much greater data storage flexibility. If you have purchased this option, please refer to the Configurable Disk Drive Supplement for descriptions of DEFAULT STORAGE operations.

Norland strongly recommends that all scan data be stored on Removable Media. This choice will minimize the amount of data that could be lost should the hard disk fail.

To access this selection screen, select **Setup** and click on the **Default Storage** command. The DEFAULT STORAGE LOCATION screen presents a choice of where future scan data is to be stored. Click on the **Fixed** button to store future scan data on the Fixed Media drive, or on the **Removable** button to store future scan data on the Removable Media volume diskettes. When satisfied with your choice, click on the **Accept Changes** button or click on the **Do Not Accept** button to abandon the operation. Until changed again, all subsequent scan data storage will routinely be on the selected disk drive.


Figure 44 - Default Storage Location Selection screen

5.4.4 PRINTER

When a printed report is requested, you are given a choice of printing immediately, or added to a printout queue. Click on the **Print Immediately** button or **Add to Queue** button. The contents of the report are assembled. If immediate printing was requested, the computer engages in the task until it is completed. If the report is to be added to the queue, the report is placed in a 'first-in-first-out' printing queue. When previously requested reports in the queue have been printed, this newly generated report is printed. To cancel the printout request, click on the **Cancel Printout** button. The initiation of printing from the print queue is made by selecting **File** and clicking on the **Begin Printouts** command.

As seen in Figure 45, Printer Setup includes all the procedures available for presetting hardcopy printing operations except type of report and choice of image rendition. The XR-26 program prints images in the same rendition as that chosen at the Analysis Results screens. The

type of report printed is the same as that seen in the Results screen. To access the other Printer setup choices, select **Setup** and click on the **Printer** command.


Figure 45 - Printer Setup screen

Select Number of Copies

We have found that many XR-26 sites prefer more than one copy of a report. XR-26 permits you to choose of the number of reports that are printed for each click the Print Report command. To select number of copies, select **Setup** and click on the **Printer** command. Click on the button for your desired number of copies. The currently preset value is highlighted and remains set until changed.

A note on the printout 'queue'. The XR-26 performs its printing on an equal basis with other XR-26 operations, so once initiated, the computer is dedicated to completing all printouts in the queue. The XR-26 program keeps track of up to 10 printout reports, provided there is sufficient space on the hard disk to hold the print files. Should a message appear saying there is no room for further reports, check to be sure the printer is turned on and is ON LINE. It then becomes necessary to wait for 'queued' reports to be completely printed before initiating other XR-26 operations.

Select Paper Length

The fan-fold printer paper used with dot-matrix printers supplied with the XR-26 system may be purchased in perforated lengths of either 11" (27.9 cm) or 12" (30.5 cm). If 11" is the XR-26 program's selected printer paper length, the printout is adjusted to be vertically centered on an 11" piece of paper. The result is that the next requested printout is centered on the next piece of fan-fold paper.

To choose the desired paper length, select **Setup** and click on the **Printer** command. Paper length choices of 11 inches or 12 inches are provided, with the currently preset paper length highlighted. Click on the desired paper length button. This selection remains set until changed.

Edit Printout Header

XR-26 printout headers may be customized by the operator as desired. The header is one line of text, up to 75 single-width characters, which is automatically centered from left to right at the top of the printed report. This line will also appear at the top of Results screens.

To change the currently stored Printout Header, select **Setup** and click on the **Printer** command. Any characters or symbols found on the computer keyboard may be entered as header information. Click on the existing header line and the editing cursor appears at the right end of the existing header. Use Dialog Edit Techniques (1.4.4) to enter or change the header information.

Leaving this screen by clicking on the **Accept Changes** button causes all of the displayed printer choices to be saved. Previously saved settings are over-written unless you click on the **Do Not Accept** button.

5.4.5 MOUSE

This Setup menu selection permits you to customize XR-26 mouse operations. As seen in Figure 46, this screen permits your choice of 'drop-down' or 'pull-down' menu operation, and your choice of mouse button double-click sensitivity. To select or make changes to these modes, select **Setup** and click on the **Mouse** command.


Figure 46 - Mouse Selection screen

Drop Down Menus

When the **Touch** button is clicked on, moving the pointer on top of the words in the Main Menu causes sub-menus to 'drop down' under the menu bar title. This is the mode as shipped from the factory.

However, beginning XR-26 instrument operators may be annoyed by seeing the sub-menus appear inadvertently. If the **Click** button is chosen, the operator must click on the Main Menu bar title before the sub-menu will appear. This choice is referred to as a 'pull-down' menu mode. As you become more familiar with XR-26 and mouse operation, the drop-down menu mode may be a more convenient choice.

Mouse Control

Routine XR-26 operations require that the mouse button be occasionally 'double-clicked' in order to rapidly progress through the program. The XR-26 program may be customized to accept your choice of double-

click response time. The average response time ('3' choice) is the mode as shipped from the factory.

If you are an experienced mouse user, you may wish to select a quicker response time, so click on the 4 or the **Fast** buttons to reduce the response time. Beginning mouse users may find that a choice of 2 or **Slow** may be easier to use. Once selected, your choice remains in effect until you change it.

Leaving this screen by clicking on the **Accept Changes** button causes the selected modes to be saved. The previous setting is over-written unless you click on the **Do Not Accept** button.

5.4.6 SCANNER SERIAL NUMBER

The Scanner Serial Number screen (Figure 47) permits you to enter the scanner serial number. This number will probably be entered by your Norland Installer. To enter the Scanner Serial Number, select **Setup** and click on the **Scanner Serial Number** command.


Figure 47 - Scanner Serial Number data entry screen

The serial number of the scanner as it exists in the computer's memory is shown. To (re)set the serial number, use the keyboard numeral keys to enter the number found on the scanner's Identification Label. When satisfied with the entry, click on the **Accept Changes** button. To discard your entries, click on the **Do Not Accept** button. Once entered, this value remains in effect until changed.

5.4.7 HIGH DENSITY DETECTION

The High Density Detection Mode is an XR-26 feature which permits the operator to select whether or not very high density measured points (those greater than about 4 g/cm² for all scan types except Whole Body BMC, which increases this limit to about 8 g/cm²) are included in the Analysis calculations. This process functions at a more primitive level than the Include / Exclude operation.

Large high density artifacts or prosthetic devices will be excluded from analysis consideration automatically if the function is enabled. However, some small metal objects may not be excluded and some small areas of bone may be excluded. High Density Detection may also be enabled or disabled at individual Modify Regions screens.

To choose a High Density Detection mode, select **Setup** and click on the **High Density Detection** command. A screen similar to Figure 48 appears.

A click on the **Exclude** button causes any scan point whose density is very high to be discarded. However, certain research protocols dictate the acceptance and analysis of a wider-than-normal range of density values. If the **Normal** button is clicked on, all points, regardless of density, are included in the Analysis calculations.


Figure 48 - High Density Detection Mode screen

Leaving this screen by clicking on the **Accept Changes** button causes the selected mode to be saved. The previous setting is over-written unless you click on the **Do Not Accept** button.

5.4.8 OTHER SETUP SELECTIONS

The remaining commands found on the Setup menu, such as X-ray Source and the Reference commands, will not be discussed here. All information concerning Reference may be found in the Reference Supplement to this Operator's Guide. The choices made in the Xray Source screen must be set only by the Norland Installer. These choices should only be changed upon direct advice from Norland Factory personnel.

WARNING: Any data acquired when the Xray Source choices are improperly set is invalid and must not be used.

5.5 TRENDING

The XR-26 provides two methods for patient comparison. Trending and Reference present a patient's bone mineral density measurement in ways which allow easy comparison to several indicators. Please note that **only** the clinician performs the diagnosis. The XR-26 supplies certain data which may assist the clinician.

One patient comparison method is called Reference, which matches the patient with sets of Reference data. Please refer to the Reference Supplement to this Operator's Guide for further information about Reference. The second method is called Trending, as described here.

Instruments such as the XR-26 measure physical bone mineral density. This property of bone is especially useful when used with other clinical parameters in assessing the status of a patient's skeleton. There is a large body of literature which suggests ways of interpreting bone mineral density measurements^{1,2}.

Comparison with Previous Measurements of Same Patient:

Probably the best interpretation of bone mineral measurements can be made by examining trends in Bone Mineral Density. If previous measurements have been made of the same patient, the current measurement may be compared with past measurements to determine whether the patient is losing or gaining bone mineral and at what rate.

¹Ross, P. D., Davis, J. W., Vogel, J. M., and Wasnich, R. D.: "A Critical Review of Bone Mass and the Risk of Fractures in Osteoporosis." Calcified Tissue International. 46:149-161, 1990.

²Melton III, J. L., Eddy, D. M., Johnson Jr., C. C.: "Screening for Osteoporosis." Annals of Internal Medicine. 112:516-528, 1990.

The rate of change of bone mineral gives the direction in which the patient's bone condition is moving, with no dependence on Reference Sets. In effect, the patient is her or his own norm.

Common sense rules were used to create the Trending and Reference operations. For example, Trending graphs will only be shown: when there is more than one scan for the patient with a matching region; there are no Reference Set(s) enabled for the region; or if "Reference/Trending" or "Trending" is the selected comparison mode. If the patient has a birthdate entered, the X-Axis label on the graph will be "AGE"; otherwise the label is "DATE". Additional information about these rules may be found in Section 7.

Before discussing Trending operations, we present several display screens and describe some of the information seen. These descriptions should assist in understanding the following detailed examples.

Figure 49 depicts a Summary Results display showing the second of four AP Spine scans for an example patient. Since the Trending operation is automatic, a trend chart is included. The values for BMD are plotted with black outlined circles. The displayed scan value is plotted in a filled black circle and may not be the most recent scan. In this example, the scan performed at age 50 is the scan being highlighted.



Figure 49 - Trend chart on Summary Results screen

To the right of the Trend Chart is a table showing the percent change in BMC that occurred since the previous scan. The short term percent change value (-17.0%) is the same as the long term value since the first and previous measurements provide identical trend results. The short term and long term percent change per year value (-16.3%) is adjusted to indicate a yearly change. See Quantities Calculated, Section 5.5.4 for further information.

The page 2 detailed results screen (Figure 50) shows the same AP Spine scan as above. Four separate charts are displayed with the tabulated measurement data.

Figure 50 - Trend charts seen on Page 2 Detailed Results screen

As you can see in Figure 51 (an enlarged portion of Figure 50), values for BMD for L2, L3, L4, and the L2-L4 column have been plotted. Again, the value for the displayed scan is a filled black circle. Since the patient history file included a birthdate, the legend at the bottom of each chart is 'AGE'.


Figure 51 - Inset of Figure 50 Trend Charts

Below each Trend Chart is a table showing similar calculated Trend information as seen on the Summary Results screen (Figure 49).

5.5.1 TREND EXAMPLE

This section describes a fictitious patient named Matilda Smith undergoing an hypothetical treatment routine. This example may provide an insight into the operation of XR-26 Trend procedures. Since all of these example screen images were obtained after a series of four measurements, each display presents all values obtained for the subject.

Patient History - Baseline Scan

Matilda Smith, age 49, presented on January 16, 1985, for scans of the AP Spine and other areas of her body in anticipation of entering a treatment program of steroid use combined with a period of inactivity.

Figure 52 is a Results display for an AP Spine scan. This first scan for the patient is intended to provide a baseline for evaluating the effects of the hypothetical treatment program in the patient's L2-L4 Spine segment.


Figure 52 - Trend charts on AP Spine Page 2 Detailed Results screen

The tabulated results values shown to the left of the Trend charts apply to the data points plotted with filled circles. Since these filled circles are the first plotted points, no previous scans have been made. Therefore, notice that no short or long term trend calculations may be made or presented.


Figure 53 - Inset of Figure 52 Trend Charts

There are several areas of interest found on this display (and printouts). Each chart 'tracks' a specific region of interest. For example, the upper left chart shown in Figure 53 (an enlarged portion of Figure 52) plots the BMD values of the L2-L4 segment against the patient's AGE. The same AGE values are used for the 'Y-axis' scale on the other Trend charts.

First Annual Scan

Following the first year of steroid treatment with inactivity, the first anniversary scan of the patient showed substantial losses of the bone mineral at sites of the spine (17%), cortical-trabecular femur (16%), trabecular femur (21%), and body (13-16%).


Figure 54 - Trend charts AP Spine Page 2 Detailed Results screen

Figure 54 shows a Results screen for the AP Spine scan taken February 1, 1986. Seen are tabulated values which apply to the data points plotted with filled circles. These filled circles are the second plotted points in the Trend charts for L2, L3, and L4 regions, plus the L2-L4 column. BMD in each region is plotted against AGE.

Following the second evaluation, the steroid treatment program was eased, a calcium supplementation program was started, and the patient was asked to take part in a twice-weekly exercise program.

Second Annual Scan

An evaluation of the third patient scan (second annual scan) showed bone mineral content values in the spine, femur, and whole body remaining essentially without change. Trend data, seen in Figure 55, indicates slight change from the previous scan data.


Figure 55 - Trend Chart on AP Spine Page 2 Detailed Results screen

Steroid treatments were discontinued after this visit, with the calcium supplementation and exercise program being continued.

Third (Final) Annual Scan

The final visit scan measurements revealed slight but significant improvements of the bone mineral in the spine and tibia, and major improvements in all sites of the femur and the bone compartment of the body.


Figure 56 - AP Spine scan Page 2 Detailed Results screen

Figure 56 trend data shows that short term %BMD changes in the L2-L4 column

(2.2) is the same as the percent per year BMD value since the previous measurement was made 12 months earlier. Note also that the long term percent change per year value has improved to -4.8.

This is the last scan series in this hypothetical study. Were the program to be continued, the Trend operations would chart up to 20 scans. Beyond 20 scans, only the 19 most recent scans plus the displayed scan would be plotted.

5.5.2 HOW TO USE TRENDING

Trending is largely automatic. If more than one scan containing similar scan type and region type exists for the patient, then a Trend chart (or charts) is always produced. The program automatically adjusts the trend chart axes to optimally plot the measurement values. This section explains the command sequences needed to obtain a Trend comparison for the selected patient.

When either Results screen is requested, the XR-26 program will present the patient's measurements using whatever information-gathering choices that were last selected by the operator. To obtain a Trending display for the current Analysis Results screen, select *Results* from the Analysis Menu title bar and click on the *Trending* command.

The display will be re-written, including one trend chart for the Summary Results screen, and several charts for the Detailed Results Page 2 screen. If the patient has a birthdate entered, the X-Axis label on the graph will be "AGE"; otherwise the label is "DATE". Appropriate Short and Long Term calculations are made and presented.

Once selected, this comparison mode remains in effect until it is changed. This means that you may select that a "Trend" chart be included in all Summary Results screens. However, these modes may be changed at any time.

5.5.3 DEFAULT CHARTED REGIONS

The XR-26 program uses the selection criteria shown in Table V to produce Results screens and printouts. The Trending information will always appear when more than one scan containing the specified values exists for the patient.

<u>Procedure</u>	<u>Value(s) Charted</u>
AP Spine	L2-L4 vs. BMD L2 vs. BMD L3 vs. BMD L4 vs. BMD Special Regions vs. BMD
Left/Right Hip	Neck vs. BMD Trochanter vs. BMD Wards Triangle vs. BMD Special Regions vs. BMD

Table V - Default Charted Regions

NOTE: Default regions are displayed from left to right and from top to bottom.

The trend information "% change" and "% change per year" will be presented only if the current scan is not the first one for the patient containing the monitored region of interest and the patient has more than one region plotted on a graph with a date of scan. The short term value indicates the percent change between the current scan and the previous scan. The long term value indicates the percent change between the current scan and the first scan shown on the graph.

If the patient has more than 20 scans of the same type, only 20 scans are plotted. If the current scan is among the most recent 20 scans of the same type for the patient, then these 20 scans are plotted.

Otherwise, the 19 most recent of that type of scan are plotted, along with the scan currently being analyzed.

5.5.4 QUANTITIES CALCULATED

When two or more serial measurements exist for the selected patient, a table of short and long term trend values is presented. The formulae used to derive these values are:

$$\text{Short} = \frac{(\text{current value} - \text{previous value}) \times (100.0)}{\text{previous value}}$$

$$\text{Short/yr} = \frac{(\text{short}) \times (365.25)}{\# \text{ days between current and previous measurement}}$$

$$\text{Long} = \frac{(\text{current value} - \text{first value}) \times (100.0)}{\text{first value}}$$

$$\text{Long/yr} = \frac{(\text{long}) \times (365.25)}{\# \text{ days between current and first measurement}}$$

This completes the description of Trending operations for the XR-26 system.

5.6 PRECISION AND ACCURACY

This section of the Operator's Guide discusses the interpretation Norland places on the terms "precision" and "accuracy", as well as our use of the precision "CV" specifications.

PRECISION

Precision is the degree to which the same value is obtained when a measurement is repeated. In-vivo precision in the XR-26 is the degree to which the XR-26 gives the same bone mineral value when a measurement is repeated at the same site on the same subject.

Precision is affected by a number of factors, some of which are inherent in the instrument, some of which have to do with the in-vivo nature of the measurement, and some of which are due to operator technique. For example, the physical phenomenon known as counting statistics will produce some imprecision in measurements of a phantom using the best of operator techniques.

Added to that for in-vivo measurements is the imprecision introduced by the fact that the subject may move slightly during the scan (if only in breathing) and that there may be some extra-skeletal calcification which changes from one measurement to the next. Finally, the operator is required to use careful judgement in selecting the bone region of interest for analysis, and this human judgement introduces some further imprecision.

Precision is normally specified as standard deviation (SD), a well-accepted statistical quantity, or as coefficient of variation (CV), which is simply the standard deviation expressed as a percent of the measured quantity.

Precision specifications given in this guide and in the XR-26 results displays and printouts are for the typical situation of (1) a properly operating instrument, (2) a healthy subject with good bone mineral content who is endeavoring to be still during the scan, and (3) an experienced operator using good analysis techniques, as described in this guide.

Imprecision may be reduced to any desired degree by averaging a number of repeated measurements. Then the precision of the mean value is improved by a factor of the square root of the number of measurements. This technique is often used, especially on phantoms, to eliminate the effects of imprecision when evaluating instrument accuracy.

ACCURACY

Accuracy is the degree to which a measurement value estimates the actual value of the quantity being measured.

Although the accuracy of a single measurement includes the imprecision of that measurement, it is customary to specify instrumental accuracy with the imprecision effect removed, as discussed above. Sometimes this is called "accuracy of the mean," because taking the mean or average of many measurements effectively removes the imprecision. The accuracy specifications for the XR-26 found in this guide are stated in terms of "accuracy of the mean".

The accuracy of measurements in instruments such as the XR-26 depends primarily on the accuracy of the calibration of the instrument. Bone mineral measurement instruments have traditionally been calibrated with reference to an "ash study", a means of determining the amount of bone mineral in an excised bone which has been scanned. There are sources of error in ash studies which have led to wide differences in calibration in instruments from different manufacturers.

In order to minimize or eliminate this inter-manufacturer variation in the recent x-ray-based instruments, several instrument manufacturers have agreed to adopt a phantom constructed of calcium hydroxyapatite and epoxy as an industry standard for bone mineral measurements. The calibration of the Norland XR-26 is based on this standard. Each XR-26 system (including its own calibration standard) is calibrated at the factory to be within 1.0% of the accepted value of the standards cited above.

When successive scans are performed on a patient using different XR-26 instruments, any patient measurement evaluation based on apparent bone mineral density changes must include the potential error due to the instrument's accuracy specification. This potential error should be treated as described in the USE OF THE PRECISION (CV) SPECIFICATIONS, later in this section.

5.6.1 PRECISION AND ACCURACY FACTORS

Precision and accuracy of the XR-26 measurements may be adversely affected by the following factors:

1. Patient Variability

The precision estimates given by the instrument are for healthy normal subjects. These estimates are for subjects approximately 20 cm thick and with spinal bone density of about 1 g/cm².

Poorer precision can be expected for subjects who:

- Are thicker than 20 cm, or
- Have less bone mineral than 1 g/cm².

2. Operator Technique

It is essential that every operator of the XR-26 follow the procedures in this guide for:

- Positioning the patient and marking the scan area,
- Keeping the patient still during the measurement, and
- Properly and repeatably defining the regions of interest.

Failure to do so can negatively affect intraoperator and interoperator precision.

3. External Radiation

Operation of the XR-26 in close proximity to external radiation sources may affect measurement results and should be avoided. External radiation sources include:

- Other x-ray generating devices,
- Devices and procedures involving radionuclides,
- Radionuclide storage areas,
- Patients containing radionuclides.

In general, if the background radiation level remains below 1 millirad per hour during the scanning procedure, the XR-26 will function properly.

5.6.2 USE OF THE PRECISION (CV) SPECIFICATIONS

An instrument such as the XR-26 estimates bone mineral content with a certain degree of uncertainty or imprecision. The clinician should take this into account when interpreting the values obtained. For example, when a

patient is re-measured after a period of time, a change in the bone mineral values may be due to an actual change in bone mineral or it may be just a variation due to imprecision.

The coefficient of variation (CV) specifications given in this guide and included in the XR-26 results displays and printouts give an estimate of the amount of imprecision that may be expected in the measurements. In order to be confident, at the 95% level, that a change in bone mineral has actually occurred, the measured change must be at least

$$2\sqrt{2} * \left(\frac{CV}{100}\right) * VALU$$

where VALU is the measured value and CV is the coefficient of variation given for that value¹.

As an example, an annual A-P Spine measurement of a subject yields a measured BMD value of 0.803 g/cm². The C.V. value for L2-L4 A-P Spine measurements is 1.0. To have 95% confidence that an actual change in bone mineral occurred since the previous measurement, the measured change for the subject must be at least 0.0227 g/cm².

A discussion of imprecision and clinical decision-making is given in a paper by Dr. Robert Heaney of Creighton University².

¹ after LeBlanc, A.D., et. al.: "Precision of Dual Photon Measurements." Journal of Nuclear Medicine. 27:1362, 1986.

² Heaney, Robert P.: "En recherche de la difference (P<0.5)." Bone and Mineral. 1:99-114, 1986.

5.7 BONE MEASUREMENT RESULTS

5.7.1 QUANTITIES ESTIMATED

The quantities estimated by measurements made with the XR-26 are defined as follows:

Bone Mineral Content (BMC)	The total mass of bone mineral within the specified region of interest, expressed in grams (g).
Bone Area (AREA)	The two-dimensional area of bone within the specified region of interest, expressed in square centimeters (cm ²).
Bone Mineral Density (BMD)	The areal density of bone within the specified region of interest; i.e., the bone mineral mass divided by the bone area; expressed in grams per square centimeter (g/cm ²).
Bone Length (LGTH)	In regions of interest which consist of a "slice" of an essentially cylindrical bone (e.g., femoral neck, spinal column), the dimension of the slice along the axis of the bone; expressed in centimeters (cm).

5.8 MECHANICAL AND ELECTRICAL HARDWARE

Norland recommends that the scanner be powered-off when not in use for periods longer than two hours. This will lengthen the service life of the x-ray source components. Instrument warm-up is not required before the performance of calibration or patient scans.

5.8.1 U. S. REGULATORY LABELS

The following illustrations portray the U.S. mandated warning and information labels which are affixed to the XR-26 instrument. All labels are shown actual size and will be found in the locations shown on Drawing Number 388D289D.

Laser Class Warning Label (378D210)

**DANGER
LASER
RADIATION
WHEN OPEN
AVOID
DIRECT
EYE
EXPOSURE**

Non-interlocked
Panel Label
(378D582)

**LASER
APERTURE**

Laser
Aperture
Label
(378D581)

**WARNING: THIS X-RAY UNIT MAY BE
DANGEROUS TO PATIENT AND OPERATOR
UNLESS SAFE EXPOSURE FACTORS AND
OPERATING INSTRUCTIONS ARE OBSERVED.**

X-Ray Warning Label (388D151)

NORLAND W6340 HACKBARTH ROAD
FT. ATKINSON, WI. 53538
MODEL XR-26 X-RAY BONE DENSITOMETER
XR-26 BEAM LIMITER
MANUFACTURED/TESTED []
MODEL 388A350 S/N []
THIS UNIT CONFORMS TO ALL APPLICABLE
STANDARDS UNDER SUBCHAPTER J,
RADIOLOGICAL HEALTH; RELATIVE TO THE
FEDERAL FOOD, DRUG AND COSMETIC ACT.

Beam Limiter Cert. Label
(388D255)

NORLAND W6340 HACKBARTH ROAD
FT. ATKINSON, WI. 53538
MODEL XR-26 X-RAY BONE DENSITOMETER
XR-26 SCANNER
MANUFACTURED/TESTED []
MODEL 388A038 S/N []
THIS UNIT CONFORMS TO ALL APPLICABLE
STANDARDS UNDER SUBCHAPTER J,
RADIOLOGICAL HEALTH; RELATIVE TO THE
FEDERAL FOOD, DRUG AND COSMETIC ACT.


Scanner Cert.
Label (388D256)

NORLAND W6340 HACKBARTH ROAD
FT. ATKINSON, WI. 53538
MODEL XR-26 X-RAY BONE DENSITOMETER
X-RAY SOURCE
MANUFACTURED/TESTED []
(REF. SOURCE LABEL FOR MODEL & S/N)
THIS UNIT CONFORMS TO ALL APPLICABLE
STANDARDS UNDER SUBCHAPTER J,
RADIOLOGICAL HEALTH; RELATIVE TO THE
FEDERAL FOOD, DRUG AND COSMETIC ACT.

Source Cert. Label (388D263)

NORLAND W6340 HACKBARTH ROAD
FT. ATKINSON, WI. 53538
MODEL XR-26 X-RAY BONE DENSITOMETER
XR-26 HOST COMPUTER
MANUFACTURED/TESTED []
(REF. COMPUTER LABEL FOR MODEL & S/N)
THIS UNIT CONFORMS TO ALL APPLICABLE
STANDARDS UNDER SUBCHAPTER J,
RADIOLOGICAL HEALTH; RELATIVE TO THE
FEDERAL FOOD, DRUG AND COSMETIC ACT.

Host Computer
Cert. Label
(388D264)


Laser Manual Shutoff (388D165)


Laser Manual Shutoff (388D395)


Other labels affixed to the XR-26 scanner:

'Hospital Grade' label

Fuse Replacement label

High Voltage Caution label

LINE VOLTAGE SELECTION	FUSE NORTH AMERICA 1/4 X 1 1/4 250V (3AG)	FUSE EXCEPT NORTH AMERICA 5X20 MM 250V
100/120 V~	6.25A SLOW	6.3 A (T)
220/240 V~	4.0A SLOW	4.0 A (T)
FREQ: 50/60 Hz~		INPUT POWER: 600 VA MAX.

Fuse Rating label

MANUAL LASER SHUTOFF LABEL

LABEL - 388D395 - (END OF SHUTOFF ARMS)
OR
LABEL - 388D165 - (END OF SHUTTER KNOB)

NON-INTER LOCKED PANEL LABEL
378D582
(UNDER ARM COVER)

LASER APERTURE LABEL
378D581 (UNDER SIDE OF
LOWER COVER)

LASER LABEL
378D210
(COLUMN FRONT)

XR-SERIES HOST COMPUTER
CERTIFICATION LABEL
388D264
(REAR OF MODEL 502, 70,
60, 555X PS/2, OR
NORLAND COMPUTER)

XR-SERIES HOST COMPUTER
CERTIFICATION LABEL
388D264
(FRONT OF MODEL 60 PS/2)

'INVESTIGATIONAL DEVICE' LABEL -
388D424 OR 'ANIMAL USE' LABEL -
388D423 LOCATIONS (3)

BEAM LIMITER
CERTIFICATION LABELS (2)
388D255, 388D444, OR 387D053
(OUTSIDE OF COVER; UNDER COVER, SIDE OF FILTER CASTING)

XR-SERIES SOURCE CERTIFICATION LABEL
388D263
(SIDE OF SOURCE)

FOR EXPORT ONLY LABEL
388D351 (FOR FOREIGN
SHIPMENTS ONLY)
(APPROXIMATE LOCATION)

X-RAY WARNING LABEL
388D151
(END OF UNIT)

SCANNER
CERTIFICATION LABEL
387D052

This drawing is the confidential property of NORLAND CORPORATION. Any reproduction of this material or its dissemination is prohibited unless specifically authorized in writing.

6345	H	2-28-91	Multiple Label IN changes	REV 01	
6463	A	1-11-91	388D264 Replaced by 387D052	REV 02	
6359	F	8-19-90	Two New Cert Labels	REV 03	
6344	A	8-7-90	Man II ECN	REV 04	
6182	A	3-9-90		REV 05	
6108	A	2-9-90	ADDED 388D264 OR 555X	REV 06	
4000	B	1-25-90	Added 388D151 Label To Draw	REV 07	
-	A	1-7-79	RELEASED	REV 08	

NORLAND CORPORATION <small>A DIVISION OF CORNING CORPORATION NORLAND DRIVE, FORT ABERDEEN, NY 14850 TEL: 716-833-4250 FAX: 716-833-4251</small>	
REGULATORY LABEL LOCATIONS PROJECT: XR-26	
REV: 8	DATE: 10-2-89
BY: [Signature]	CHK: N.T.S.
APP: [Signature]	DATE: 11/1/89
REV: 1	DATE: 11/1/89
REV: 1	DATE: 11/1/89

5.8.2 SCANNER DRIVE SYSTEM

The scanner incorporates a precision X-Y drive capable of resolving position within 0.1 mm in each direction, and of repositioning within 0.2 mm. Arm movement speed is variable from 1.0 mm/sec to 80.0 mm/sec.

Motion in each direction is provided by a stepper motor and a timing belt drive. The stepper motors are operated by specially designed hardware which controls speed and distance of motion and keeps track of position. The two directions are completely independent, meaning that the scanner may move in any angular direction at any speed.

The stepper control hardware is controlled by the scanner microprocessor, which in turn receives commands from the IBM computer.

The communication between the control/analysis work station and the scanner is handled via an opto-isolated three-wire RS-232C interface. XR-26 software expects that the host computer is dedicated to the program. Use of other software programs on the computer may cause problems with XR-26 operations. Prior to installing other application programs on the host computer, please contact your local Norland representative to discuss the matter.

5.8.3 PHOTON DETECTION SYSTEM

The detection system consists of two NaI (TI) crystal scintillation detectors coupled to proprietary amplification and counting circuitry.

5.8.4 SCANNER CONTROL PANEL

The Scanner incorporates a touch-sensitive Scanner Control Panel. The scanner control panel is located at the outward end of the scanner arm and is depicted in Figure 57.


Figure 57 - XR-26 Scanner Control Panel

Buttons on the control panel permit the operator to turn the laser on and off and to move the scanner arm in four directions to virtually any position on the table, all while standing at the side of the patient. These functions assist in identifying certain patient anatomical locations. When the control/analysis unit acquires these specific locations, it computes the proper scan area, then carries out the scan.

Also located on the scanner control panel are indicator lights and an Emergency OFF key. The **EMERGENCY OFF** key is used to instantly stop all scanner operations and turn off the x-ray beam and laser. Separate indicators illuminate whenever the X-RAY source has

power applied to it, whenever the SHUTTER is open, and whenever the LASER is on.

5.8.5 COMPUTER SYSTEM

The computer used with the XR-26 instrument is a member of the popular IBM⁴ Personal Computer family. XR-26 Software is currently configured for use with the IBM-PS/2 personal computer family (Models 50Z, 55SX, 60, 70, or 80), with minimum 1 Mbyte RAM memory, minimum 30 Mbyte Hard Disk, 1.44 Mbyte floppy disk drive, MS-DOS V3.3 or 4.0, 80287 or 80387 Math co-processors, VGA color graphics, and high resolution display.

The XR-26 Instrument is equipped with a high-resolution graphics printer and a computer work station.

⁴ Trademarks and Registered Trademarks owned by International Business Machines, Inc.

5.9 PREVENTIVE MAINTENANCE

The XR-26 is largely maintenance-free. Status of electronic modules is checked at every start up or when the self-test procedures built into the computer software are initiated. However, there are mechanical parts of the system which do need some periodic maintenance in order to remain in top working order.

5.9.1 SCANNER

Do the following at least once each month:

Clean guide rails: The guide rails are the round steel rods inside the scanner which guide the motions of the scanner arm and source carriage assembly. There are two sets of rails, one set for each direction.

To gain access to the rails, turn off and unplug the scanner from the power source. Gently lift the front edge of the patient pad/table top and detach it from the scanner. Wipe the guide rails clean with a dry, lint-free cloth or paper towel. The arm/carriage assembly may have to be moved about in order to clean the entire length of the rails. When the rails are clean, replace the patient pad and reconnect the scanner power cord.

Clean patient pad and positioning blocks: For general cleaning, an ordinary vacuum cleaner or 'DustBuster' equipped with a stiff-bristled sweep attachment serves the purpose. To remove lint,

hair, or small particles, a short-nap or adhesive-covered lint-brush is suitable.

For stubborn dirt, gently scour soiled places on the patient pad or positioning blocks with a stiff-bristled brush. The use of a foaming-type upholstery cleaner approved for use with flame-retardant fabric is acceptable. Blot up excess cleaner with a cloth or paper towel. Avoid letting the cleaning solution drip off the patient pad into the scanner base. Avoid getting the cleaning solution into the zippers of the positioning blocks. Avoid using products which may promote static build-up.

It is recommended that ordinary hospital sheets and pillowcases be used to preserve the cleanliness of the patient-contact surfaces. Should laboratory animals be the subject of investigation, it is prudent to cover the scan surface with a thin moisture barrier. Materials such as rubber sheets or polyethylene no thicker than 5 mil will not interfere with the scan results.

NOTE: There is NO maintenance required for the laser positioning system or the x-ray source. The Quality Control Procedures outlined in Section 2.3.1 (especially the daily calibration and phantom scan) insure proper operation and identification of problems, should they arise.

5.9.2 COMPUTER SYSTEM

The computer requires very little maintenance. Check the user reference manuals packed with the computer for maintenance details about the monitor, mouse, and computer.

Norland highly recommends cleaning the heads in the floppy disk drive at least once each month. To do so, obtain a head cleaning kit in a

size appropriate for the computer's disk drive from any local computer supply store. Follow the directions supplied with this kit.

It may also prove necessary to clean the CRT display screen periodically. Use an ordinary glass cleaner and lint-free towels for this task. Avoid letting the cleaner seep down inside the bezel of the monitor or into the computer.

5.9.3 HOW TO RELOCATE THE INSTRUMENT

Precautions must be taken when moving the computer. To protect it's hard disk drive, please follow the directions for relocating the computer given in the operation manual. Note that the warranty does not apply to misuse of the instrument. Do not move the computer when the power is on or if either disk drive is spinning. Damage and/or loss of data may result.

The scanner must be moved **ONLY** by qualified Norland Service personnel. Special tools and programs must be used to assure proper alignment and leveling of the scanner.

5.9.4 DIAGNOSTICS AND SERVICE

The highest quality customer assistance is provided directly from your Norland Service representative. In the event of an equipment failure or erratic measurement results, please be prepared to discuss the problem on the telephone. Telephone numbers and other information helpful for contacting Norland may be found in the Service statement at the beginning of this guide.

In general, make note of any unusual instrument operation as it occurs. We suggest that the failed operation be tried a second time. Also note

whatever conditions which may have contributed to the instrument failure.

Error Recovery

In very rare situations, it is possible for the XR-26 software or scanner to 'lock-up' or to display Error Messages. The recommended recovery procedure is to:

- A. Make a written note about what may have caused the problem, especially error codes and error messages, if any.
- B. Press the **Enter** key on the keyboard. Pressing this key causes the computer to reset itself back to DOS. Any partially-processed scan or analysis data will be lost and the XR-26 program will have to be re-started.
- C. If a problem may have been caused by the scanner, it is a good idea to turn the scanner off, then on again after a brief pause. Doing so will re-initialize the scanner's hardware and may eliminate further problems.

The least-recommended error recovery strategies are the 'cold-boot' and 'warm-boot' processes. These methods result in the absolute loss of any in-process scan or analysis data. If all else fails, to perform a 'warm-boot', press and hold both the '**Ctrl**' and '**Alt**' keys, then tap the '**Del**' key. The computer will re-boot back to the '**C>**' DOS prompt. To perform the 'cold-boot', simply turn Off the computer, scanner, and printer via their power switches, then turn them back on. XR-26 software may then be re-started for normal operations.

Fuse Replacement

All components of the XR-26 Instrument are protected by line (mains) fuses. In the event of a fuse failure in the computer or monitor, refer to the operation manual that was supplied with the unit.

The XR-26 scanner is equipped with a power switch that includes the line voltage input connector and the fuse holder. Instructions for changing the fuse is given in the Fuse Changing Instruction Sheet (388D152) which is attached to the spare fuses supplied with the XR-26 Scanner. Please refer to the label affixed near the fuse holder to determine the proper rating for the replacement fuse.

The Norland Installer configured the Scanner to operate properly on the line voltage and frequency furnished at the installation site. If the Scanner is to be operated at a different line voltage or frequency, please contact the factory for assistance in converting and moving the unit. Improper conversion is not covered by the warranty.

Supplies

Supplies for the XR-26 system may be obtained from local sources or directly from Norland (Note: Minimum order is US\$75.00). Norland part numbers for various supplies may be obtained from Norland Customer Service. Supplies may also be purchased from most office supply or computer stores, using the following guidelines:

For the IBM-PS/2 computers, purchase 3 1/2", Double Sided, High Density, (Verbatim 86269 or equivalent) diskettes.

Printer paper should be 9-1/2 inches wide and 11 or 12 inches long, fan fold, tractor feed, with perforations for bursting to 8-1/2 inch wide sheets.

Printer ribbons are available from computer and office supply stores.
They may also be ordered directly from the printer manufacturer.
An address is given in the printer manual.

This completes the descriptions of subjects related to the basic operation of the Norland XR-26 and XR-26 Mark II instruments. Supplements to this Base Operator's Guide describe all the other optional features of these bone densitometry systems.

This page reserved.